ImaCor Zura



1000 1011





ImaCor Inc 839 Stewart Avenue Suite #3 Garden City, New York 11530

Table of Contents

New in Zura 2.1.6	б
Enhanced hTEE measurement package	6
Reduced recovery period after cineloop acquisition	6
Enhanced cine comment window	6
Enhanced hTEE tutorials	6
hTEE reports	6
Expanded data transfer module (DTM) capabilities	6
Still-frame export capability	6
	-
Chapter 1: Introduction	7
User Guide Conventions	7
	/ 7
	, 7
Numbered Steps	/ 7
Fourmerated Lists	/ 7
	/ 7
	/
Crassial Terres and Assessment	/
Special Terms and Acronyms	/
Chapter 2: The ImaCor Zura Imaging System	8
Overview	8
Indications for Use	8
Contraindications	8
Transesophageal Echocardiography (TEE) Use in Critical Care Settings	8
Principles of Operation	8
Episodic Assessment	8
Non-Imaging Mode10	0
Cumulative Imaging Time10	0
About the System	0
Ultrasound Machine1	1
Ultrasound Probe	2
Adaptaflex [™] Technology12	2
Chapter 3: Getting Started	3
Prenaring the System 1	3
Connecting and Bemoving the ImaCor Disposable Probe	3
Connecting the ECG	5
Turning the System On/Off	6
System Software Controls	6
General Description of Licer Controls	7
Operator Consolo	, 7
Using the Trackhall	, Q
Monitor 1	o o
	0
Chapter A. Inc. Cox Illtracound Inc. ring Software	•
Chapter 4: Imacor Oltrasound Imaging Software	9
Overview	9
Online Help System	9
Additional Help Resources	9
System Initialization	9
Patient Information	1
Tab Order 23	3
Patient List	3
Patient List Sort Order24	4
Patient Information Screen Button Functions	4
Patient Information Tasks2	5
Add a New Patient2	5
AutoSelect a patient record26	6

Review an Existing Patient Record and Associated Cineloops	27
Import Patient Data	27
Export Patient Data	28
Querying a Modality Worklist (MWL) Server	29
Synchronizing with a DICOM Archiving Server	29
Delete a Patient	30
Preparing to Image	31
Troubleshooting	31
The Imaging Environment	
ClariTEE Plug and Play	31
Single View Mode	32
Real-Time Imaging	36
Poal-Time Imaging	36
Auto-Q Feature	
Filter Control	38
Contrast Control	
Depth Control	39
Real-Time Image Quality Controls for Color Flow Imaging	40
CF Gain Control	
CF Filter Control	
The ECG Waveform	
FCG Cursor	
FCG Markers	41
Changing the ECG Load	лл. л1
Changing the ECG Sween Speed	۰۰۰۰۰ ۹ ۱ ۸۱
Left Ventricle Trans-Gastric Short-Axis View (TGSAV)	
Depth Markers	
Freeze/Unfreeze	43
• ·	
Steering	43
Steering Image Acquisition	43 44
Steering Image Acquisition Recovery Mode	43 44 46
Steering Image Acquisition Recovery Mode Loading a Cineloop	43 44 46 46
Steering Image Acquisition Recovery Mode Loading a Cineloop Reviewing a Cineloop	43 44 46 46 47
Steering Image Acquisition Recovery Mode Loading a Cineloop Reviewing a Cineloop Cineloop Image Quality Control	43 44 46 46 47 47
Steering Image Acquisition Recovery Mode Loading a Cineloop Reviewing a Cineloop Cineloop Image Quality Control Brightness Control	43 44 46 46 47 47 47
Steering Image Acquisition Recovery Mode Loading a Cineloop Reviewing a Cineloop Cineloop Image Quality Control Brightness Control Depth Brightness Compensation (DBC)	43 44 46 46 47 47 47 47
Steering Image Acquisition Recovery Mode Loading a Cineloop Reviewing a Cineloop Cineloop Image Quality Control Brightness Control Depth Brightness Compensation (DBC)	43 44 46 46 47 47 47 47 47
Steering Image Acquisition Recovery Mode Loading a Cineloop Reviewing a Cineloop Cineloop Image Quality Control Brightness Control Depth Brightness Compensation (DBC) Contrast Control	43 44 46 46 47 47 47 47 47
Steering	43 44 46 46 47 47 47 47 47 47 47
Steering Image Acquisition Recovery Mode Loading a Cineloop Reviewing a Cineloop Cineloop Image Quality Control Brightness Control Depth Brightness Compensation (DBC) Contrast Control Depth Control Cineloop Control Buttons	43 44 46 46 47 47 47 47 48 48
Steering Image Acquisition Recovery Mode Loading a Cineloop Reviewing a Cineloop Cineloop Image Quality Control Brightness Control Depth Brightness Compensation (DBC) Contrast Control Depth Control Cineloop Control Buttons Scrolling Mode	43 44 46 46 47 47 47 47 48 48 48
Steering Image Acquisition Recovery Mode Loading a Cineloop Reviewing a Cineloop Cineloop Image Quality Control Brightness Control Depth Brightness Compensation (DBC) Contrast Control Depth Control Cineloop Control Buttons Scrolling Mode Color Flow Indicator Icon	43 44 46 46 47 47 47 47 47 48 48 48 49 49
Steering Image Acquisition Recovery Mode Loading a Cineloop Reviewing a Cineloop Cineloop Image Quality Control Brightness Control Depth Brightness Compensation (DBC) Contrast Control Depth Control Cineloop Control Buttons Scrolling Mode Color Flow Indicator Icon Exporting a Cineloop as a Movie	43 44 46 47 47 47 47 47 48 48 48 49 49 50
Steering Image Acquisition	43 44 46 46 47 47 47 47 48 48 48 49 49 50 50
Steering	43 44 46 46 47 47 47 47 48 48 48 49 50 50
Steering Image Acquisition	43 44 46 46 47 47 47 47 47 48 48 48 49 50 50 51
Steering Image Acquisition Recovery Mode Loading a Cineloop Reviewing a Cineloop Cineloop Image Quality Control Brightness Control Depth Brightness Compensation (DBC) Contrast Control Depth Control Cineloop Control Buttons Scrolling Mode Color Flow Indicator Icon Exporting a Cineloop as a Movie Saving the Current Frame as a BMP Image Measurement Tools The Tracing Tool Creating a Trace	43 44 46 46 47 47 47 47 48 48 48 49 50 50 51
Steering Image Acquisition Recovery Mode Loading a Cineloop Reviewing a Cineloop Cineloop Image Quality Control Brightness Control Depth Brightness Compensation (DBC) Contrast Control Depth Control Cineloop Control Buttons Scrolling Mode Color Flow Indicator Icon Exporting a Cineloop as a Movie Saving the Current Frame as a BMP Image Measurement Tools The Tracing Tool Creating a Trace Editing a Trace	43 44 46 46 47 47 47 47 47 48 48 48 49 50 50 51 51 53
Steering Image Acquisition	43 44 46 46 47 47 47 47 47 48 48 48 48 49 50 50 51 51 53 54
Steering Image Acquisition	43 44 46 46 47 47 47 47 47 48 48 48 48 49 50 50 51 51 53 54 54
Steering . Image Acquisition . Recovery Mode . Loading a Cineloop . Reviewing a Cineloop . Cineloop Image Quality Control . Brightness Control . Depth Brightness Compensation (DBC) . Contrast Control . Depth Control . Cineloop Control Buttons . Scrolling Mode . Color Flow Indicator Icon . Exporting a Cineloop as a Movie . Saving the Current Frame as a BMP Image . Measurement Tools . The Tracing Tool . Creating a Trace . Editing a Trace . Deleting a Trace . Deleting a Trace . The Measuring Tape Tool . Taking a Measurement	43 44 46 46 47 47 47 47 47 48 48 48 48 49 50 50 51 51 51 53 54 54 54
Steering . Image Acquisition . Recovery Mode . Loading a Cineloop . Reviewing a Cineloop . Cineloop Image Quality Control . Brightness Control . Depth Brightness Compensation (DBC) . Contrast Control . Depth Control . Cineloop Control Buttons . Scrolling Mode . Color Flow Indicator Icon . Exporting a Cineloop as a Movie . Saving the Current Frame as a BMP Image . Measurement Tools. The Tracing Tool . Creating a Trace . Editing a Trace . Deleting a Trace . Deleting a Trace . The Measuring Tape Tool . Taking a Measurement . Editing a Measuremen	43 44 46 46 47 47 47 47 47 48 48 48 48 49 50 50 50 51 51 51 54 54 55 56
Steering	43 44 46 46 47 47 47 47 47 48 48 48 48 49 50 50 50 51 51 51 53 54 54 55 56
Steering	$\begin{array}{cccccccccccccccccccccccccccccccccccc$
Steering . Image Acquisition . Recovery Mode . Loading a Cineloop . Reviewing a Cineloop . Cineloop Image Quality Control . Brightness Control . Depth Brightness Compensation (DBC) . Contrast Control . Depth Control . Cineloop Control Buttons . Scrolling Mode . Color Flow Indicator Icon . Exporting a Cineloop as a Movie . Saving the Current Frame as a BMP Image . Measurement Tools . The Tracing Tool . Creating a Trace . Deleting a Trace . Deleting a Trace . Deleting a Trace . Deleting a Measurement . Editing a Measurement . Deleting a Measurement . Viewing Measurement .	$\begin{array}{cccccccccccccccccccccccccccccccccccc$
Steering . Image Acquisition . Recovery Mode . Loading a Cineloop . Reviewing a Cineloop . Cineloop Image Quality Control . Brightness Control . Depth Brightness Compensation (DBC) . Contrast Control . Depth Control . Cineloop Control Buttons . Scrolling Mode . Color Flow Indicator Icon . Exporting a Cineloop as a Movie . Saving the Current Frame as a BMP Image . Measurement Tools. The Tracing Tool . Creating a Trace . Deleting a Trace . Deleting a Trace . Deleting a Trace . Deleting a Measurement . Editing a Measurement . Deleting a M	$\begin{array}{cccccccccccccccccccccccccccccccccccc$
Steering . Image Acquisition . Recovery Mode . Loading a Cineloop . Reviewing a Cineloop . Cineloop Image Quality Control . Brightness Control . Depth Brightness Compensation (DBC) . Contrast Control . Depth Control . Cineloop Control Buttons . Scrolling Mode . Color Flow Indicator Icon . Exporting a Cineloop as a Movie . Saving the Current Frame as a BMP Image . Measurement Tools . The Tracing Tool . Creating a Trace . Editing a Trace . Deleting a Trace . Editing a Trace . Deleting a Measurement . Deleting a M	$\begin{array}{cccccccccccccccccccccccccccccccccccc$
Steering . Image Acquisition . Recovery Mode . Loading a Cineloop . Reviewing a Cineloop . Cineloop Image Quality Control . Brightness Control . Depth Brightness Compensation (DBC) . Contrast Control . Depth Control . Cineloop Control Buttons . Scrolling Mode . Color Flow Indicator Icon . Exporting a Cineloop as a Movie . Saving the Current Frame as a BMP Image . Measurement Tools . The Tracing Tool . Creating a Trace . Editing a Trace . Deleting a Trace . Deleting a Trace . Deleting a Measurement . Editing a Measurement . Deleting a Measurement . Deleting a Measurement . Viewing Measurements . Viewing Measurements on Cineloops . Fractional Area Change (FAC) Calculations . The FAC Wizard .	$\begin{array}{cccccccccccccccccccccccccccccccccccc$

Trace Overlay Mode	62
Superimposed Comparison	02
	02
Split View Mode	02
Loading Two Cincloops in Split View Mode	64
Viewing a Cincleon and Peal Time Images in Split View Mode	
Suppression of Dayback	05
Synchronized Playback.	00
Activating and Deactivating Synchronized Playback	6/
Playback Speed Controls	6/
	69
The Color Box	
B-Mode Priority	/2
Setting b-mode priority	/3
Velocities	73
Color Maps	74
CF Overlay in Playback Mode	75
Printing	77
Single View Mode	77
Split View Mode	77
Configuring the System	79
General	79
Imaging Defaults	79
Measuring and Tracing Settings	79
Institution	80
Reset Button	80
Acquisition	80
Date-Time	81
DICOM	81
Language	83
Reports	84
hTEE Reports	84
Version	86
Software Update	86
Download the Software Upgrade and Request a License File	87
Upgrade the Zura Software	87
Help	87
Chapter 5: Safety	91
Disposable Single-Use ClariTEE	
The ALARA Principle	92
Applying ALARA	92
ALARA Operator Control	92
	92
Mechanical Index	92
ALARA Automatic Controls	92
Twenty Minute Software Interlock	22
Watchdog Monitor	92 02
Transoconbaggal Echo Imaging (TEE) Potential Complications	92
Electrical Safety	
Electromagnetic Compatibility (EMC)	95
Piocompatibility Safety	
Acoustic Output	
Global Maximum Devated ISDTA, ML and Ti Values	95 0F
A coustic Moosurement Provision and Uncertainty	
Labeling Symbols	96
Chapter 6: ImaCor Approved Devices	97
Chapter 7: Troubleshooting and Maintenance	98

Troubleshooting	
Maintenance	
Zura System Cleaning Recommendations	
LCD Monitor Cleaning	
Operator Console Cleaning	
Handle and Umbilical	100
Other Recommended Maintenance	100
Power Supply Cord	100
Handle and Umbilical	
Chapter 8: Master Message Listing - Alphabetical	101
Chapter 9: References	108
Chapter 9: References	108
Chapter 9: References. Measurements (Area and Distance). Sources of Measurement Errors	108 108 108
Chapter 9: References. Measurements (Area and Distance). Sources of Measurement Errors Accuracy	108 108108108
Chapter 9: References. Measurements (Area and Distance). Sources of Measurement Errors Accuracy Methodology of Assessing Accuracy.	108 108108108108108108
Chapter 9: References. Measurements (Area and Distance). Sources of Measurement Errors Accuracy Methodology of Assessing Accuracy. Chapter 10: Specifications.	
Chapter 9: References. Measurements (Area and Distance) Sources of Measurement Errors Accuracy Methodology of Assessing Accuracy Chapter 10: Specifications. Chapter 11: Glossary.	
Chapter 9: References. Measurements (Area and Distance). Sources of Measurement Errors Accuracy Methodology of Assessing Accuracy. Chapter 10: Specifications. Chapter 11: Glossary Terms	



Obelis s.a.

Bd. Général Wahis 53, 1030 Brussels Belgium +32 2 732 59 54



ImaCor Inc. 839 Stewart Ave Suite 3 Garden City, NY 11530 Phone: 516-393-0970 Fax : 516-393-0969

DECLARATION OF CONFORMITY

(Annex II.3 of the Directive 93/42/EEC)

ImaCor Inc. declares conformance with the relevant provisions of the aformentioned EC directive. This declaration is issued under the sole responsibility of ImaCor Inc.

Manufacturer :]
	8

ImaCor Inc 839 Stewart Ave Suite 3 Garden City, NY 11530 USA

Product Identification : ImaCor TEE Imaging System consisting of

Ultrasound Machine Zura, Zura-EVO, EVZ100, ZT-1000 with handle model UMB-1000

ClariTEE probe

<u>Classification :</u> The Ultrasound Zura, Zura-EVO machine model is classified as Class IIa in accordance with rule 10 of Annex IX. The ClariTEE probe is classified as Class IIa in accordance with rule 5 of Annex IX

European Uni	on Representative :	Obelis s.a. Bd. General Wahis 53, 1030 Brussels, Belgium Tel : +32(0)2 7325954		
European Notified Body :		TÜV Rheinland LGA Products GmbH Tilllystrasee 2 D-90431 Nürnberg Tel.: +49 221 806 2004 Fax: +49 228 806 3935 Identification number: 0197		
Quality Manag	<u>ement System :</u> Ac	cording to EN ISO 13485 2003 + AC		
Signed	Richard Lanzillotto	ned by Richard Lanzillotto ward Lanzillotto, e=ImaCor illottosimacorinczon, ceUS Date		

Richard Lanzillotto Director of Regulatory Affairs

New in Zura 2.1.6

Enhanced hTEE measurement package

Zura 2.1.6 enables you to perform three types of hTEE measurements:

- FAC (LVEDA, LVESA)
- RV/LV ratio (RVEDA, LVEDA)
- SVC collapsibility index (SVC max, SVC min)

All three measurements can be taken directly. Fractional area change (FAC calculations can also be performed with the assistance of the FAC Wizard, which guides you through the process of completing FAC calculations.

Reduced recovery period after cineloop acquisition

Recovery time has been reduced from 9 seconds to 3.3 seconds.

Enhanced cine comment window

View type (TGSAV, ME4CH, SVC, or other) can be recorded in the cine comment window. The selected view type will be displayed in the cineloop load dialog window and will be used by the system to automatically display the correct measurement package (FAC RV/LV, SVC). See page 45.

Enhanced hTEE tutorials

hTEE tutorials now include video clips.

hTEE reports

Zura now includes hTEE reporting capabilities, including pre- and post-assessment reporting and a summary report. See page 84.

Expanded data transfer module (DTM) capabilities

USB external drives can now be formatted as DTMs, accommodating large import/export operations.

Still-frame export capability

Cineloop frames can be saved as BMP images on a USB thumbdrive).

Chapter 1: Introduction

This user guide covers the set-up and operation of the ImaCor Zura system, which is intended for use only by qualified clinicians. It does not provide instruction in sonography, cardiology, or other clinical practices. Please carefully read this guide before using the ImaCor Zura[™] imaging system. Pay special attention to the information contained the Safety section, beginning on page 91, as well as text marked with the **A Caution** symbol.

This user guide covers the set-up and operation of the ImaCor Zura system, which is intended for use only by qualified clinicians. It does not provide instruction in sonography, cardiology, or other clinical practices.

User Guide Conventions

lcons

Caution: Describes a procedure or precaution necessary to prevent injury to the patient or damage to the system.

Safety Feature: Highlights a safety feature.

The system a feature unique to the ImaCor Zura system

Text Conventions

NUMBERED STEPS

Numbered steps are formatted as shown below. The steps **must be performed** in the order listed. *Example–*

There are three steps in the process:

- 1 Step one
- 2 Step two
- 3 Step three

ENUMERATED LISTS

Enumerated lists are formatted as shown below. In this case, the numbers indicate quantity, not order.

Example-

- The dropdown menu lists three options:
- 1. Option one
- 2. Option two
- 3. Option three

BULLETED LISTS

Bulleted lists are formatted as shown below. They represent a group of associated items, but do not imply quantity or order.

Example-

- System features include:
- Feature
- Feature
- Feature

TERMS

Button:An onscreen software controlKey:A key found on the ECG console keyboardKnob:A knob found on the ECG console keyboard

Special Terms and Acronyms

A glossary of terms appearing in this guide begins on page 111. A list of acronyms and definitions begins on page 111.

Chapter 2: The ImaCor Zura Imaging System

Prescription use: For use by qualified clinicians only.

TimaCor Innovation

Miniaturization of the TEE probe permits assessment of cardiac function episodically over an extended period of time, as the smaller probe should not require the use of general anaesthesia.

FimaCor Innovation

The ImaCor single-use disposable probe can remain in place for up to 72 hours, enabling episodic assessment of the patient's cardiac function over that period at intervals determined by the physician.

Important: The Zura system is not a continuous monitoring system.

Overview

- The ImaCor miniaturized TEE probe (ClariTEE[™]) enables direct visualization of cardiac size and function, and is designed specifically for episodic assessment in the critical care environment.
- Miniaturization of the probe permits an extended maximum dwell time of 72 hours.
- The single-use probe is provided sterile and is disposable.

Indications for Use

The ImaCor system is intended for the **episodic assessment** of cardiac function using transesophageal echocardiography (TEE). It is indicated for use in clinical settings, including long-term settings such as the ICU, for an indwelling time period not to exceed 72 hours.

Contraindications

The ImaCor Zura system is contraindicated as follows:

- It is not for pediatric use (patients less than 18 years of age).
- **Warning:** Patients with known issues such as esophageal or stomach varices, obstructive esophageal pathology, recent surgery, or radiation therapy in the esophageal or gastric area, should be evaluated by a physician prior to having a TEE procedure.

Transesophageal Echocardiography (TEE) Use in Critical Care Settings

The ImaCor Zura system addresses in the critical care setting the need for a miniaturized transesophageal echocardiography (TEE) probe in assessing important cardiac parameters that influence hemodynamics. The new ultrasound probe (ImaCor ClariTEE[™]) provides direct visualization of cardiac size and function, allowing intensive care clinicians to conduct episodic assessments of cardiac preload and left ventricular ("LV") systolic function over an extended period.

Principles of Operation

Miniaturization of the ClariTEE probe enables direct visualization of the left ventricle episodically over time. Direct visualization of the left ventricle permits improved assessment of intravascular volume status and cardiac function.

The ClariTEE disposable probe can remain in place for up to 72 hours, enabling episodic assessment of the patient's cardiac function over that period at intervals determined by the physician. Episodic assessments typically occur up to six times over a 24-hour period.

Episodic Assessment

The ImaCor Zura system is not a continuous monitoring device. It is intended to be used to conduct episodic assessments of the patient's cardiac function. Normal use in this setting should consist of a maximum of six episodic assessments over a 24-hour period with maximum intubation time not to exceed 72 hours. The duration of a typical assessment will be in the range of 5–15 minutes. Hence, the typical cumulative imaging time for each patient is in the range of 1.5 hours to 4.5 hours during a 72-hour indwelling time period.

Table 1 Episodic Assessment Overview

Anticipated episodic assessment frequency	Six episodic assessments per 24 hours
Typical episodic assessment time	5–15 minutes
Maximum probe dwell time	72 hours
Typical cumulative imaging time during a 72-hour indwelling time period	1.5–4.5 hours
Maximum cumulative imaging time during 72-hour time period	6 hours



Fig. 1 Episodic assessment process flow

Safety Feature

During non-imaging intervals, no energy is delivered to the patient and the probe remains in a relaxed, unflexed position. As a result, the mucosal tissue is not subject to any mechanical or thermal stress.

A Caution

The Zura system should be in imaging mode only when an episodic assessment is being conducted.

Non-Imaging Mode

The ClariTEE probe is designed for a maximum dwell time of 72 hours. Episodic assessments typically account for a fraction of total dwell time; therefore, during the majority of probe dwell time no imaging is taking place. During non-imaging intervals, no energy is delivered to the patient and the probe remains in a relaxed, unflexed position.

Cumulative Imaging Time

Cumulative imaging time is the total time over the course of probe dwell duration that the Zura system is in imaging mode for the purpose of episodic assessment. The ClariTEE probe delivers energy to the patient only while in the imaging mode. Also, only in imaging mode is the ClariTEE probe mechanically flexed to obtain the tissue contact required to capture an image.

About the System

The ImaCor Zura system consists of four main components:

- 1. Ultrasound machine
- 2. ClariTEE probe
- 3. Ultrasound handle
- 4. ImaCor ultrasound imaging software



Fig. 2 ImaCor Zura Ultrasound Machine

Ultrasound Machine

The ImaCor Zura system includes an ultrasound machine optimized for use with the ImaCor ClariTEE probe (Figure 2). The machine contains a liquid crystal display (LCD), an operator console, and all the required system firmware and hardware except for the disposable ClariTEE probe.

A Caution

The ImaCor ClariTEE is for use **only** with the ImaCor Zura System

Flex the probe only when imaging.

Disconnect the probe from the handle when not conducting an episodic assessment.

The ClariTEE probe is designed for single-use only.

The maximum dwell time for the ClariTEE probe is 72 hours.

Ultrasound Probe

The ClariTEE is a miniaturized disposable single-use probe optimized for extended dwell-time, making it ideal for use in longer-term clinical settings such as the ICU. During episodic assessments, the probe tip is flexed upward.

The indwelling portion of the ClariTEE probe is detachable from the control handle. The detachable handle enables one machine to serve multiple patients.

ADAPTAFLEX™ TECHNOLOGY

The ClariTEE probe features ImaCor AdaptaFlex technology. This technology changes the flex point of the probe's adaptive section to accommodate the patient's anatomy. AdaptaFlex technology enables the operator to obtain an optimal view across patients of varying sizes.



Fig. 3 ImaCor ClariTEE[™] probe, umbilical, and handle

Chapter 3: Getting Started

Preparing the System

Plug the system into a hospital grade or equivalent receptacle (outlet).

Connecting and Removing the ImaCor Disposable Probe

- 1 Locate the dedicated matching connector (ZIF connector) for the locking umbilical cable connector.
- **2** Push the umbilical cable connector into place (Figure 4). The umbilical cable connector is large; you will need one hand to hold it and the other to turn the locking lever.



Fig. 4 Push the umbilical connector

3 When ready to attach the probe, remove it from its sterile packing. Remove the cover from the disposable probe connector as shown (Figure 5).



Fig. 5 Remove the cover from the disposable probe

4 Position the probe connector and handle as shown in Figure 6.



Fig. 6 Position the probe connector and handle in hand

5 Align the probe with the handle.



Fig. 7 Align probe and handle

6 Insert the blue connector tab into the distal end of the handle (Figure 8).



Fig. 8 Insert connector tab into distal end of the handle

7 Gently press the probe connector over the handle; press the tab on the umbilical connector to assist with engagement (Figure 9). An audible click will indicate the probe is properly connected.



Fig. 9 Press the probe connector over the handle

8 With the disposable probe locked into place, check articulation by moving level in each direction. The distal tip of the probe should flex in either direction.
 Note: Expect flex of 90 degrees in the forward direction (anteflex) and 20 degrees minimum in reverse direction (retroflex).



Fig. 10 Check articulation of the probe

Connecting the ECG

A three-lead ECG cable is provided for optional connection of the machine ECG to the patient. The ECG provided with the Zura is not for diagnostic purposes and is provided only to assist the detection of images at end systole and end diastole. After connecting the patient electrodes, connect the ECG cable to the machine using the circular twist-lock receptacle located next to the printer.



Fig. 11 Connect the ECG lead

Turning the System On/Off

 Turn on the power supply by pressing the rocker switch on the back of the unit to the ON position, marked I. (The OFF position is marked O.)



Fig. 12 ImaCor Zura on/off switch

2 Power up the computer by pressing the pushbutton switch directly behind the ECG connector shown in Figure 11 and sharing the same housing.



Fig. 13 Power pushbutton switch

System Software Controls

See Figure 14 on page 17 for a detailed console drawing.

Table 2 System Software Console Controls

Control	Function
Alphanumeric keys	Enable user to enter comments or data into text boxes
"C" key	Toggles Color Flow imaging mode on and off. This can also be accomplished with the on-screen Color Flow button.
View	Toggles the display view from Single View mode to split screen mode (two loops displayed side by side). For addi- tional information, see Single View mode, page 31 and Split View Mode, page 63.
Print	Prints the displayed frozen (still) ultrasound image. See Printing, page 77.
Acquire	Stores a three-second loop of a real-time ultrasound images and ECG data (if any). See Image Acquisition, page 44.
Freeze	Toggles the Freeze/Unfreeze of a real-time image. See Freeze/Unfreeze, page 43.
Trace	Selects trace mode, which enables the user manually to trace areas. See Tracing Tool, page 51.
Measure	Selects measure mode, which enables user to measure height and width distances; e.g., in the LV cavity. See Measuring Tape Tool, page 54.
Overlay	Sustains the display of traced LVEDA and LVESA contours over other frames in a loop. See Trace Overlay Mode, page 62.
Load	Loads a cineloop into the display. See Loading a Cineloop, page 46.

ImaCor Ultrasound imaging software performs seven basic functions:

- 1. Recording and updating patient information
- 2. Real-time imaging
- 3. Cineloop acquisition
- 4. Cineloop enhancement
- 5. Cineloop playback
- 6. Cineloop evaluation
- 7. FAC calculation

GAIN	 Real-time mode: controls overall gain Playback mode: controls brightness
	See Gain Control, page 36.
TGC	Real-time mode: controls TGC
	 Playback mode: controls DBC
	See Time Gain Compensation Control (TGC), page 40.
FILTER	Real-time mode: controls filter and contrast
	 Playback mode: controls contrast
	See Filter Control, page 40.
DEPTH	Real-time mode: controls displayed image depth
	 Playback mode: controls displayed image depth
	See Depth Control, page 48.
"N" key	Loads the next cineloop without going through the Load
	Loop dialog.
"P" key	Loads the previous cineloop without going through the
	Load Loop dialog.

Table 2 System Software Console Controls (Continued)

General Description of User Controls

Operator Console



Fig. 14 ImaCor Zura operator console

Using the Trackball

The trackball is a pointing device. It serves the same functions as a mouse, but accomplishes them in different ways. A mouse must be moved from one location to another to reposition the cursor or to scroll, whereas a trackball is rotated in its socket while the device in which it is housed remains stationary. Trackballs permit fast, continuous scrolling and need not be repositioned.

As a rule, trackballs are controlled with the thumb and fingertips, allowing for finer control than a mouse, which is operated with the hand and wrist.

- Roll the trackball in the direction you wish to move the cursor.
- Roll the trackball left or right to scroll

You must be in scrolling mode to scroll through images with the trackball. Left click on the image to enter scrolling mode; left click again to exit.

Like a mouse, the trackball device includes right and left buttons.

- Click the left button on the trackball to select a software button or a menu item.
- Click the right button on the trackball to view a context-sensitive menu.

Monitor

The LCD monitor contains brightness controls. These can be adjusted to the user's specifications. Other monitor controls should left at factory default settings.

Trackballs permit fast, continuous scrolling and need not be repositioned, In addition, trackballs allow for finer control than a mouse.

Chapter 4: ImaCor Ultrasound Imaging Software

Overview

ImaCor Ultrasound Imaging Software v2.1.6 performs seven basic functions:

- 1. Recording and updating patient information
- 2. Real-time imaging
- 3. Cineloop acquisition
- 4. Cineloop enhancement
- 5. Cineloop playback
- 6. Cineloop evaluation
- 7. hTEE measurement

Online Help System

The complete User Manual is available via electronically through the imaging software's searchable Help System. To access the Help system, press and hold the key, then press the H key, or right click and select "Help" from the global context menu.



Fig. 15 Sample Help screen

Additional Help Resources

Additional assistance in the form of tooltips and live technical support is also available. See page 87 for details.

System Initialization

See Troubleshooting (page 31) and Master Message Listing (page 101) for an explanation of system errors and suggested solutions. The first screen displayed is the system initialization screen. During initialization, configuration files are loaded into memory and the system software conducts a diagnostic self-test. Informational messages are displayed if a problem is detected.

See Troubleshooting (page 31) and Master Message Listing (page 101) for an explanation of possible failures and suggested solutions.

A full-power test is run after a patient examination ends abruptly; i.e., when an exam is ended without the user pressing the End Exam button. Causes for an abrupt exam termination include a power outage, a system powered down during an exam, and a software malfunction. The full-power test adds about 40 seconds to the initialization sequence.



Fig. 16 System initialization screen – normal initialization sequence

Patient Information

Basic patient data are entered using the Patient Information window (Figure 17). After starting the ImaCor Zura system, the Patient Information screen is the first interactive window displayed.

Ima ImaCo	r Lab										11/16/2011 05:10:02 PM
	Patient inf	formation:	Ê								
	guid:	{D122BDB	E-D4E4-4D6C-9E	09-0A978FF5	D954}						
	Patient ID:	1202230					Notes:				<u> </u>
	Last Name:	DOE									
	First Name:	JOHN									
	DOB:		mm/dd/yyyy								V
	Weight:	lb	oz	Imperia	al 🔤	-					
	Height:	ft	in				Physicia	n: Firs	t Name	Last Name	
	Sex:	Unknown	•				Accessio	on #:			
		Last Name	First Name	DOB	Sev	Physician	Cine	Date Created	Accession #		Enable Eiltering
	0020110	DANTON	1ENNIEER	1000	U	Filysicidii	0	11/16/2011	Accession #	{D4ACABE	
	0291102	JORDAN	STEVEN		Ŭ		16	11/3/2011		{593219BB	
	1202230	DOE	JOHN		U		• 1	11/3/2011		{D122BDBE	
	2216254	WILBURG	MICHAEL	01/01/1902	U		0 46	10/7/2011		{462EE118	
Imp	ort Exp	orta				New Patie			Next	Configure	e Shutdown

Fig. 17 Patient information screen, information editing enabled (see "Patient Information Screen Button Functions," page 24).

The fields detailed in Table 3 comprise the Patient Information screen.

Table 3 Patient Information Screen Data Fields

Field Name	Status	Function	Notes
GUID	System generated	A unique patient identifier	Read-only; i.e., the information in this field is displayed, but cannot be changed
Patient ID	Optional	Enter local patient identification num- ber, if applicable	32 characters maximum
Last Name	Mandatory	Enter patient's last name	To ensure consistency, first and last name information is always displayed and recorded in upper- case characters. 32 characters maximum

Field Name	Status	Function	Notes
First Name	Mandatory	Enter patient's first name	To ensure consistency, first and last name information is always displayed and recorded in upper- case characters. 32 characters maximum
DOB	Optional	Enter patient's date of birth	Format is mm/dd/yyyy; e.g., 06/02/1956.
Weight	Optional	Enter patient's weight	The drop-down menu to the right of the Weight field enables you to select Imperial or metric units.
			If you change units after entering the patient's weight, the mea- surement is converted.
			The option you select becomes the system default.
Height	Optional	Enter patient's height	The drop-down menu to the right of the Height field enables you to select Imperial or metric units.
			If you change units after entering the patient's height, the mea- surement is converted.
			The option you select becomes the system default.
Physician	Optional	Enter the physician's first and last name	Each field can accommodate 32 characters.
Sex	Optional	Enter patient's sex	The drop-down menu provides three choices: 1. Male 2. Female 3. Unknown
Accession #	Optional	Enter accession number (can also be imported from an MWL server)	The accession number is an iden- tifier for an Imaging Service Request. The Zura system allows a single accession number per patient. The accession number is embedded in DICOM files that are sent to the archiving server (when configured).
Notes	Optional	Enter additional comments	This freeform text field can accommodate 256 characters.
Enable filtering	Optional	Enable/disable patient list filtering	When the patient list filter is turned on (the option box is checked), the characters you type in the first and last name fields limit the names displayed in the patient list to those with matching character strings.
			For example, if you type "Sm" in the last name field and "J" in the first name field, only patients whose last name begins with the string "Sm" and first name begins with "J" will be displayed: • Smith, John • Smith, Josephine • Smithers, Joachim

Table 3 Patient Information Screen Data Fields (Continued)

Tab Order

Use the Tab key to advance through the enterable fields in the Patient Information screen. The tab order is:

- Patient ID
- Last Name
- First Name
- DOB
- Weight
- Unit System
- Height
- Sex
- Notes
- Physician First Name
- Physician Last Name
- Accession Number

Patient List

A list of patient records appears at the bottom of the Patient Information screen.

ID 🔺	Last Name	First Name	DOB	Sex	Physician	Cine	Date Created	Accession #	GUID
00000003	DOE	JANE	06/11/1948	F	JESSE, JONES	02	1/10/2008		{51882AEA
00000055	MELVILLE	ROBERT	12/08/1960	U	PAUL, DAWSON	02	1/20/2011	12345678	{D657B971
00000104	VERNON	MIKE	01/02/1955	U	JON, DELERUE	01	4/28/2011		{BF67216A
00000105	JACKSON	STANLEY	01/12/1954	U	ROB, SCHMITT	01	4/28/2011		{6260E5E1
00000117	MC QUEEN	JAMES	11/02/1940	М	JONES, JESSE	05	8/20/2010		{AD84F7D7
1									

Fig. 18 Sample patient list

The list includes the following information:

- Local identification number, if applicable
- Patient's last name
- Patient's first name
- Patient's date of birth
- Patient's sex
- Physician name
- Cine the number of cineloops recorded for the patient (maximum 100)
- Synchronization state icon (under Cine) indicates the state of synchronization with the archiving server. There are three states:
 - No sync (no cineloop for the selected patient has been sent to the DICOM server)
 - Partial sync (at least one cineloop for the selected patient has been sent to the DICOM server)
 - Complete sync (all cineloops for the selected patient have been sent to the DICOM server)
- Patient record creation date
- Accession #
- Patient GUID

PATIENT LIST SORT ORDER

You can set the sort order of the patient list by clicking on the appropriate column heading.

You can set the sort order of the patient list by clicking on the appropriate column heading. For example, click on the Last Name column heading to display the patient list alphabetically by last name. The active column is indicated with an arrow.





Arrow indicating active column

Fig. 19 Patient list column headings

Patient Information Screen Button Functions

Figure 20 shows the patient information screen buttons in their default state.



Fig. 20 Patient information screen buttons - default state

The Import and Export buttons are enabled when an ImaCor Data Transfer Module (DTM) is connected. See page 28 for information in importing and exporting patient data.

The Sync button is enabled when a DICOM Archiving Server is active (see "Archiving Server Configuration," page 82). The Sync button allows you to send cineloops as DICOM files to a DICOM server via a network connection. Right-click the Sync button to view a server synchronization log.

The new atomic DICOM synchronization feature permits you to send selected patient records to the archiving server. Previously, all patients were archived in a synchronization operation.

Server Sync Log 🛛 🛛 🗙
Last Sync Date: 05/23/2011 - 01:35:54 PM
Synchronization Errors: {AD84F7D7-7259-4CB9-B459-4EE98BBB6A2F},08-19-2010,20100819-183450.cine,1223 {AD84F7D7-7259-4CB9-B459-4EE98BBB6A2F},08-20-2010,20100820-001054.cine,1223
Number of skipped loops: 32 Rules not complied with: PatientD: PatientDDB: 0 PatientSex: 0 AccessionNum: 25 CineComment: 0 MaxNumberDaily: 0
Number of loops sent out: 1
ОК

Fig. 21 Representative server synchronization log

Active buttons are slightly dimmed to minimize visual distractions. The New Patient, Next, Configure, and Shutdown buttons are active but slightly dimmed. This is the default state of active buttons in the Imacor system; it is intended to minimize visual distractions.

Number	Button	Function
1	Import	Import patient information. Enabled only when a Data Transfer Module (DTM) is inserted in a USB port.
2	Export	Export patient information. Enabled only when a Data Transfer Module (DTM) is connected and a patient record is selected from the patient list.
3	Sync	Synchronize cineloops with a DICOM Archiving Server. See DICOM configuration.
4	New Patient	Clears the Patient Information screen, including the GUID, enabling the user to create a new patient record. The New Patient button is the only way to clear the GUID.
5	Next	Advances to the imaging preparation screen. Mandatory first and last name fields must have been completed.
6	Configure	Accesses the Configuration dialog screens; see "Configur- ing the System," page 79, for details.
7	Shutdown	Shuts down the system.

on
•

Table 5 Standard ImaCor Button States

Button	Button State
	Disabled
	Default
	Highlighted (rollover)
	Selected (trackball click)

When you roll the trackball over an active button, the button is highlighted as shown in Table 5. Clicking on an active button changes the color from blue to orange, also shown in Table 5. This sequence is standard throughout the ImaCor Zura system.

Patient Information Tasks

ADD A NEW PATIENT

- 1 If data from an existing patient are displayed in the Patient Information screen, press the New Patient button.
- 2 Enter the new patient information in the appropriate fields. Press the Tab key to advance from one field to the next; Press Shift-Tab to move backwards. The First Name and Last Name fields are mandatory; all others are optional.
- **3** When you have finished entering information, press the Next button to begin imaging preparation.

Autosave Feature

Once a single character is entered in both the First Name and the Last Name fields, the autosave feature records it and all subsequent characters in the patient database. No user action is required save information entered in the Patient Information screen.

Patient Information Validation

Information entered in the Patient Information screen is validated before you enter the imaging environment. If an error is detected, the field is highlighted.

Patient Edit Lock Feature

The Patient Edit Lock feature prevents unintentional changes to critical patient data fields. Once a patient record is created, the patient ID, last name, and first

name fields are locked. To edit these fields in an existing record, you must unlock them. The lock icon next to the Patient Information screen title displays the current state of a selected record.



The icon on the left indicates that a selected patient record is unlocked and that all fields are available for editing. The icon on the right indicates that patient records are locked and editing of the patient ID, last name, and first name fields is disabled.

The lock icon works like a toggle switch. Click the icon to lock and unlock a selected patient record. When you unlock a record, a confirmation message is displayed.

AUTOSELECT A PATIENT RECORD

Whenever a probe is connected to the Zura system, the system software checks to see if the probe has been used before. If the probe has been used in a previous patient exam, the system reads the patient data encoded in the probe—the GUID and the patient's first and last names—then searches for a match in the patient database. If a match is found, the system will, at your direction, automatically select the record and populate the data fields. The user prompt is shown below in Figure 22.



Fig. 22 Autoselect prompt

If you reject the autoselect option, the system software turns off the probedetection mechanism for the duration of the patient exam or until the probe is disconnected.

Patient Records Created on a Remote System

If patient data has been encoded in the probe, but no matching record is found in the patient database, the Zura system software automatically creates a patient record based on the encoded information. The system then displays the autoselect prompt.

This situation arises when a patient's initial examination is completed on one Zura system, but a subsequent assessment uses a different Zura. For example, a patient's initial exam may take place in an emergency room setting, while later examinations are completed in the intensive care environment. These departments likely use different Zuras, and Zuras are stand-alone systems. Information entered in one patient database does not appear in databases associated with other Zura systems. This is also true for cineloops.

A system-generated patient record will contain just the patient's GUID and first and last names. Other patient information, as well as associated cineloops, must be imported from the remote Zura system (see "Import Patient Data," page 27). Information that is corrected in or added to a system-generated patient record will not be reflected in the original record.

UPDATE A PATIENT RECORD

Use the Patient Information screen to include additional information in a patient record or to correct existing information.

- 1 When the Autoselect prompt is displayed, click the Yes button; alternatively, you may manually locate the patient's record in the patient list. If you choose the latter, click on the patient's name to select the record. The patient's personal information is displayed in the Patient Information screen data fields.
- **2** Tab to the field you wish to modify.
- 3 Enter the new information or correct the existing data.

If you access the Patient Information screen from the imaging screen during an exam, press the OK button to return to the imaging environment.

REVIEW AN EXISTING PATIENT RECORD AND ASSOCIATED CINELOOPS

You must select an existing patient record when you wish to review previously recorded cineloops or acquire new ones. A cineloop is a three-second loop of ultrasound images and ECG data.

- 1 When the Autoselect prompt is displayed, click the Yes button; alternatively, you may manually locate the patient's record in the patient list. If you choose the latter, click on the patient's name to select the record.
- 2 Press the Next button to advance to the imaging preparation environment.
- **3** For information on viewing an existing cineloop, see "Loading a Cineloop," page 46.

IMPORT PATIENT DATA

Patient information and cineloops can easily be exchanged between Zura systems. To import or export information, you must have a Data Transfer Module (DTM), which is a jump drive specially formatted for transferring ImaCor Zura system records.

Importing data into the ImaCor Zura system is a three-step process:

- 1 Insert the DTM into one of the two USB ports on the face of the console. When the DTM is connected, the Import button becomes active and the system displays a message indicating that the Export feature is enabled.
- 2 Click the Import button. The DTM Patients dialog box is displayed. Click the name of the patient whose records you wish to import. Pressing the Shift button while clicking enables you to select multiple records.

M patients	ultinle DTM pati	ents to import					
ID -	Last Name	First Name	DOB	Sex	Physician	Cine	Date Created
0000003	DOE	JANE	06/30/1934	F	JONES, J	12	10/15/2008
00000055	MELVILLE	ROBERT	05/23/1936	М	DAWSON, P	4	1/10/2008
(1
						Import	Cancel

Fig. 23 DTM patient display

3 Click the Import button. Do not disconnect the DTM while the data is being imported. Figure 24 shows the Importing progress bar.

Note: Patient data on the DTM is encrypted. A validation (checksum) of the DTM data is performed by the system when a DTM is detected.

Importing		×					
Patient:	{4655EA3C-CBF3-4350-967B-4DF81EBCE564}						
File:	20120624-081302.dcf						
Importing patient data Do not disconnect DTM.							
Cancel							

Fig. 24 Importing patient data progress bar

Duplicate Records

If you are importing a record that already exists on the target machine, a warning message is displayed, as shown in Figure 25.

- If you wish to overwrite the record on the target machine, click the Yes button. Overwritten records cannot be retrieved.
- If you wish to retain the record on the target machine, click the No button. The patient record on the DTM will not be transferred.

Conflict Detected							
	Conflict:	Last Name					
	Patient:	FRED ARDEL					
	r daene.						
	Select value below t	resolve the conflict					
	Select value below (Diesolve die connict.					
	• Local:	ABDEL					
	O DIM:	AREL					
	· • • • • •	HDEE					
		OK					

Fig. 25 Patient record overwrite warning

EXPORT PATIENT DATA

Exporting data from an ImaCor system is a three step process:

- 1 Insert the DTM into one of the two USB ports on the face of the console. The instructional message shown in Figure 26 is displayed.
- 2 Click the name of the patient whose records you wish to export. Pressing the Shift button while clicking enables you to select multiple records.
- **3** Click the Export button.

Data Transfer Module (DTM) connected	X
The Export and Import functions have been enabled.	
To export a patient to the DTM, select the patient in the patient list and click on the Export button.	
ок	

Fig. 26 DTM confirmation

Notes:

- Patient data is encrypted before it is stored on a DTM.
- If the "Export DICOM Images" option is enabled in the Configuration window, the operation will export ImaCor cine files as well as DICOM multiframe images.
- The ImaCor decryption utility must be used before DICOM images can be imported and viewed on a DICOM workstation.

QUERYING A MODALITY WORKLIST (MWL) SERVER

If an MWL server is configured and networked to a Zura system, you can query the server for a selected patient. Data elements that are already present on the Zura system are not imported from the MWL server. For example, if the date of birth for a selected patient is already recorded on the Zura system, the birthdate will not be imported.

- 1 Select a patient from the patient list.
- 2 Right-click and select "Check MWL Server...".

If the MWL server query returns at least one item, you can import the information from a selected MWL item into the local patient record. See DICOM Conformance Statement for details of imported fields.

Note: The accession number can be imported from the MWL server.

SYNCHRONIZING WITH A DICOM ARCHIVING SERVER

If a DICOM archiving server is configured and networked to a Zura system, you can send acquired cineloops to it as multi-frame DICOM files

Transmit cineloops to the server by clicking the Sync button in the patient screen. After communication with the server is verified, the Zura system sends cineloops that do not yet reside on the server and that satisfy user-defined Sync rules, if any. Examples of user-defined rules follow:

- Require cineloop comment
- Require patient ID
- Require patient DOB
- Require accession number
- Require patient sex ("unknown" not acceptable)

Note: Cineloops imported from another Zura system will not be sent to the archiving server. Cineloops can transmitted to an archiving server only by the Zura system that acquired them.

The Sync Log feature enables you to check the status of the last synchronization operation, displaying transmission errors and the count of cineloops not transmitted because of rule noncompliance. To view the Sync Log, right click the Sync button, which displays the Show Synchronization Log menu.

The new atomic DICOM synchronization feature permits you to send selected patient records to the archiving server. Previously, all patients were archived in a synchronization operation.

DELETE A PATIENT

- Locate the patient's record in the patient list.
- Right click the patient's name, then choose the Delete option as shown in Figure 27, *or*—
- Press the DEL or Backspace console key.

To delete multiple patients, select the patient names and delete as described above.

Ima ImaCoi	r Lab										12/1/2011 10:04:50 AM
	Patient in	formation:									
	guid:	{110F66	62-621E-4D98-B9D0	C-7B111E0D9A	A1}						
	Patient ID:	212423	42				Notes:				<u> </u>
	Last Name:	SMITH									
	First Name:	MICHA	EL								
	DOB:		mm/dd/yyyy								
	Weight:		b oz	Imperia	il 🔤	•					
	Height:	f f	t in				Physicia	an: Firs	st Name	Last Name	
	Sex:	Unknow	n 💌				Accessi	on #:			
	ID 🔺	Last Name	First Name	DOB	Sex	Physician	Cine	Date Created	Accession #	GUID	Enable Filtering
	123412344	WALLACE	JOHN		U	, ,	17	7/27/2011		{3C2BE258	
	21242342 00000117	SMITH MC QUEEN	MICHAEL Check Warklist Delete	11/02/1940	M	JONES, JESSE	0	11/8/2011 8/20/2010		{110F662 {AD84F7D7	
Imp	ort Exp	oort Sy	nc			New Patie	ent		Next	Configure	e Shutdown

Fig. 27 Deleting a patient record

When the confirmation message is displayed (Figure 27), press the Yes button to delete the patient record. Deleted records cannot be retrieved.



Fig. 28 Deleting a patient record - confirmation message

Preparing to Image

Once you've completed the Patient Information screen, the Imacor Zura system begins a series of self checks in preparation for imaging.

The Imaging Preparation screen has three buttons: Back, Next, and End Exam.



Fig. 29 Imaging Preparation screen buttons

3. Imaging Preparation Screen Buttons and Function

Number	Button	Function
1	Back	Returns you to the Patient Information screen.
2	Next	Advances you to the Imaging screen
3	End Exam	Halts patient exam and displays the Patient Information screen

Troubleshooting

If the ImaCor Zura system encounters a problem, an error message is displayed on the Imaging Preparation screen. See "Chapter 7: Troubleshooting and Mainte-nance," page 98 for a complete listing of error messages and their meaning.

The Imaging Environment

The imaging environment is where you will view, acquire, optimize, and compare real-time images and cineloops.

There are two basic screen modes in the imaging environment: Single View and Split View. Single View mode displays one real-time or cineloop view. See page 32 for information on Single View mode. Split View mode enables you to compare two cineloops side-by-side or one cineloop beside real-time imaging. See page 63 for information on Split View mode.

ClariTEE Plug and Play

The ClariTEE probe is detected and validated from the Imaging Screen. This enables you to disconnect/reconnect without having to exit the Imaging Screen to validate the reconnected probe.

Single View Mode

Single View mode is generally used for real-time monitoring and FAC calculations. To review and compare saved cineloops with live images, or to compare two saved cineloops, see "Split View Mode," page 63.

See Table 6 on page 33 for information on Single View mode features.



Fig. 30 Single View mode imaging screen

Number	Feature	Function	Notes
1	Probe meter	Displays the probe's unique serial number and a progress bar indicating the time remaining on the probe (0–72 hours).	
2	Patient name	Patient's first and last name.	
3	Frame title	Alerts user that real-time imaging is taking place.	Indicates current view. Dis- plays Real Time when viewing real-time images. When a cineloop is loaded, displays the cineloop title; i.e., the date and time of acquisition.
4	hTEE report	Opens the hTEE Report window	
5	hTEE buttons	Displays the hTEE Views tutorial	
6	Current date and time		
7	Tx frequency	Displays the current ultra- sound transmit frequency. The values are 6 MHz and 7 MHz for B-mode imaging and 5 MHz for Color Flow imaging. The default B-mode is transmit 7 MHz. It can be changed to 6 MHz for improved penetra- tion with a slight decrease in resolution.	To toggle between B-mode frequencies, click within the Frequency area, bounded with a blue rectangle, to display the frequency context menu. Select the desired frequency from the menu. The default B-mode fre- quency setting can be changed in the Configuration dialog box. The Color Flow frequency can- not be changed.
8	Steering feature cursor	Enables the user to look 25 degrees to the left or right without physically moving the probe.	 Steering is always available from the keyboard during real-time imaging. Press > to pan right on image Press < < to pan left on image Press < < to pan left on image Steering also available from the trackball unless measuring or tracing mode is ON. Left-click on the right of the image to pan right Left-click on the left of the image to pan left The steering angle increments in five degree steps.

Table 6 Imaging Screen Features

Table 6 I	maging Screen Featu	ires (Continued)	
Number	Feature	Function	Notes
9	Grayscale	Used in adjusting monitor brightness and contrast.	The grayscale is a vertical gra- dient bar on the imaging screen that allows the user to adjust the LCD's brightness to an optimal setting for the ambient lighting. The monitor should be adjusted to display as much of the grayscale range as possible.
			Grayscale is also used to set the Color/B-Mode Priority when imaging with CF or replaying a CF cineloop.
10	Color Map	A 256-color scale used to map velocities in the Color Flow data.	Left-click the color map to step through the five available color maps.
11	Velocity display	Low/Medium/High	Only visible in CF mode
12	Depth display	Displays the imaging depth. The default imaging depth is 15 cm.	Always ON
13	ECG waveform		
14	Heart rate	Expressed in beats per minute (bpm)	
15	Frame rate	Expressed in frames per second (fps)	Standard frame rate for B- mode imaging is 24 fps
			The B-mode acquisition frame rate is 50 fps.
			The CF-mode imaging frame rate varies width colorbox width. The lowest frame rate is 6 fps for the widest colorbox; the highest frame rate is 20 fps for the narrowest color box.
			The CF acquisition frame rate is double the standard CF frame rate with a maximum of 30 fps.
			Recovery frame rate is 6 fps. For more information on recovery mode, see page 46.
			Note: The frame rate shown represents the number of frames per second displayed on screen. This number is lower than the fps captured during acquisition due to the time required to process and render images on screen.
16	ECG lead	Select Lead I, II, or III	
17	Color indicator		Only visible during real-time CF imaging
18	ECG image		
19	Institution (hospital) name and location		



Fig. 31 Single View mode layout buttons

Table 7 Imaging Screen Buttons

Number	Feature	Function	Notes
1	Freeze/Unfreeze	Toggles between Freeze and Unfreeze mode	Freeze: Halts imaging on the last frame viewed and internally saves the last 50 frames.
			The 50-frame buffer enables you to scroll backward to view previous frames. When scroll- ing, the fps rate display indi- cates the speed at which you are reviewing frames.
			Use the trackball to scroll through the buffered frames.
			Use the playback buttons to click through the sequence one frame at a time.
			Unfreeze: Resumes real-time imaging.
			See "Freeze/Unfreeze," page 43
2	Color Flow	Toggles Color Flow imaging on and off.	During real-time imaging, tog- gles Color Flow imaging mode on and off. This can also be accomplished with the "C" console key.
3	Acquire	Starts a 3-second acquisition of real-time images and ECG.	See "Image Acquisition," page 44. You may acquire a maximum of 100 cineloops per patient.
4	Patient	Displays the Patient Informa- tion screen, enabling you to edit patient information with- out ending an exam. To return to imaging view, click the OK button.	If the Patient button is clicked during an exam, real-time imaging is automatically fro- zen although the exam is still in progress.
			After patient information is edited and the OK button pressed, you are returned to the imaging screen and the exam resumes.
			If a cineloop was being viewed, it remains loaded.
			If real-time imaging was taking place, the frozen 50-frame real-time buffer is displayed.
5	Print	Sends the current image to a local printer, along with FAC data if applicable.	See "Printing," page 77
6	Single/Split View	Toggles between Single View and Split View modes.	See "Single View Mode," page 32 and "Split View Mode," page 63
Number	Feature	Function	Notes
--------	-----------	---	--
7	FAC Calc	Starts the FAC Wizard, which guides the user through the steps necessary to complete an FAC calculation.	See "The FAC Wizard," page 58
8	Measure	Displays Measure pop-up menu, allowing user to select Area, Distance, or Off option.	
9	Load	Loads a cineloop.	See "Loading a Cineloop," page 46
10	Configure	Enables user to customize several aspects of the system configuration, including: • Application • Date-Time • DICOM • Language • Version • Help	"Configuring the System," page 79
11	End Exam	Terminates patient exam and returns the user to the patient information screen.	A patient exam is completed by clicking the End Exam but- ton or by selecting a different name from the patient list.

Real-Time Imaging

Real-Time Image Quality Controls for B-Mode Imaging

The real-time image quality controls for B-mode enhance frames that are being captured by the ultrasound hardware. The real-time control parameters are automatically saved when a patient is being imaged and are restored when the patient is re-examined.

Default values can be restored for each control by pressing the Heart console key (\bigcirc) in combination with the appropriate image quality knob. For example, pressing \bigcirc plus the Gain knob restores the default values for Gain and CF Gain.

Note: Default values can be customized in the Configuration window; choose the General tab.

GAIN CONTROL

The Gain Control knob, located on the console, regulates the overall gain of the ultrasound system. Gain controls the overall brightness of an image. It is used in combination with TGC, as described below.

Press or rotate the Gain Control knob to display the Image Quality Control display with the gain controls active (Figure 32).

The gain range is 0-20; the default setting is 10. The current setting is displayed on the right side of the gain gauge.



Fig. 32 Gain Control display

TIME GAIN COMPENSATION (TGC)

TGC control is an index in a family of preset TGC curves. It simplifies the usual TGC slider controls.

Since ultrasound waves are attenuated, or absorbed, as they pass through tissue, the waves reflected from distant areas are weaker than those from the areas near the transducer. Without TGC, the image would have a light-to-dark gradient from the near field to the far field. TGC enables you to compensate for the loss of ultrasound waves in the far field by amplifying their signal. The result is even brightness across the entire field of view.

Press or rotate the TGC control knob to view the Image Quality Control display with the TGC controls active. The TGC range is 0–19 (Figure 33); the TGC1 default setting is 5. TGC2 control is similar but controls the rolloff on the TGC curve. This is useful to attenuate the brightness in the far field. Press the TGC knob to toggle between TGC and TGC2. The TGC2 range is 0–16.



Fig. 33 TGC controls display

AUTO-Q FEATURE

Auto-Q is a "smart" Gain+TGC function that automatically adjusts the Gain, TGC1, and TGC2 settings until an optimal image is obtained.

To activate the Auto-Q feature, press the "Q" key on the console or left-click the on-screen Auto-Q button during live imaging.



Fig. 34 Auto-Q feature

Note: Image quality is subjective; you may prefer Gain and TGC settings other than those selected by Auto-Q.

FILTER CONTROL

The Filter control is a despeckling (noise reduction) filter that enhances the contrast of the preprocessed ultrasound images. The Filter control range is 0–9, with 0 specifying no filtering and is adjusted in one-step increments.

Press or rotate the Filter knob to display the Image Quality Control display with the Filter controls active (Figure 35).



Fig. 35 Filter control display

CONTRAST CONTROL

If the Filter control is displayed, press the knob once to toggle to Contrast control (Figure 36). If another control is displayed (e.g., Gain, TGC, or Depth) or if the Control display is not active, press the Filter knob twice to access Contrast control. The settings are from 0 to 10 in one-step increments.



Fig. 36 Contrast control display

DEPTH CONTROL

The Depth control determines the depth in centimeters of the displayed image.

Press or rotate the Depth knob to display the Image Quality Control display with the Depth controls active (Figure 37). The Depth control range is 6–15 cm and is adjusted in 1 cm increments.

The default imaging depth is 15 cm.

The Zura system automatically saves the last imaging depth for the selected patient and restores that depth when the patient is again selected for imaging.



Fig. 37 Depth control display

True Depth feature

Images captured at depths less than 12 cm are magnified to fill the same display area as those obtained at depths \ge 12 cm.



Fig. 38 True Depth feature—the same image at depths of 7 cm, 10 cm, and 12 cm (left to right).

Real-Time Image Quality Controls for Color Flow Imaging

Color flow imaging quality is managed by two real-time image quality controls: CF Gain and CF Filter. B-mode controls continue to be available in CF mode.

CF GAIN CONTROL

Press or rotate the Gain Control knob on the console to access CF Gain. Press the Gain knob to toggle between B-mode gain and CF gain. The CF gain range is 0–45.



Fig. 39 CF gain control

CF FILTER CONTROL

The CF Filter control range is 1–21 and is adjusted in one-step increments. The Filter knob toggles between CF filter, B-mode filter, and B-mode contrast.

Press or rotate the Filter knob to display the Image Quality Control display with the CF Filter controls active (Figure 35).



Fig. 40 CF filter control

The ECG Waveform

The ImaCor Zura system enables you to customize the ECG waveform display.



Fig. 41 Example ECG waveform with markers

ECG CURSOR

The ECG cursor is a gray vertical bar (playhead) that indicates the current frame in the waveform being viewed, as shown in Figure 41. When reviewing frames, the ECG cursor is highlighted as it passes over an ECG marker, as shown in Figure 42.



Fig. 42 ECG waveform with highlighted cursor

ECG MARKERS

By default, an ECG marker is automatically placed at the peak of the R wave (see Figure 41). Marker placement is usually reliable, but can be inaccurate if:

- there is no ECG signal
- the signal is too low for R-wave detection
- the R-wave detection software misidentifies the R wave, which can occur with extremely atypical ECG waveforms

In these instances, ECG markers can be entered and edited manually.

Adding an ECG Marker

- 1 Move the ECG cursor to the desired position.
- **2** Right click the trackball button.
- 3 Choose the Add ECG Marker option from the pop-up menu (see Figure 46)

Deleting an ECG Marker

- 1 Move the ECG cursor until it is highlighted over the marker you wish to remove.
- 2 Right click the trackball button.
- **3** Choose the Delete ECG Marker option from the pop-up menu (see Figure 46)

If you delete all ECG markers on a cineloop, the ImaCor System software will attempt to detect the R wave peaks and place ECG markers accordingly the next time the cineloop is loaded.

CHANGING THE ECG LEAD

- 1 Right click the trackball button.
- 2 Choose the desired ECG Lead from the pop-up menu (see Figure 44)

CHANGING THE ECG SWEEP SPEED

- 1 Right click the trackball button.
- 2 When the pop-up menu opens, choose from the two sweep-speed options: 25 mm/s and 50 mm/s, as shown in Figure 46.

INVERTING THE ECG WAVEFORM

- 1 Right click the trackball button.
- 2 Choose the Invert ECG option from the pop-up menu (see Figure 46).

Figure 45 below shows an inverted ECG waveform.



Fig. 43 Global menu



Fig. 44 Change ECG lead option



Fig. 45 Inverted ECG waveform with markers (markers larger than actual size)



Fig. 46 ECG menu options

Left Ventricle Trans-Gastric Short-Axis View (TGSAV)

In Single View mode, real-time images of the left ventricle in the TGSAV are displayed in the center of the screen (see page 45). In Split View mode, images are centered in the right and left halves of the screen, as shown on page 64. The Zura system is optimized for the TGSAV, but other views associated with TEE imaging may be obtained.

In both Single and Split Views, images with a depth of 12 cm are life size; i.e., a distance of 1 cm on screen represents an actual distance of 1 cm.

- The new True Depth feature magnifies images with depths from 6 cm-11 cm to fill the same life-size area as a 12 cm image.
- Images with depths from 13 cm-15 cm are scaled down.
- The maximum image depth is 15 cm.



Fig. 47 Example left ventricle image in the TGSAV

DEPTH MARKERS

To display depth markers on the ultrasound image:



Fig. 48 Display depth marker option

- 1 Right click the trackball button.
- 2 Choose the Depth Markers option from the global context menu (see Figure 48). The interval between two adjacent depth markers is 1 cm.

Depth markers are superimposed on the ultrasound image on the right of the image, as shown in Figure 49. Depth marker setting is persistent across exams and system reboots.



Fig. 49 Real-time image with depth markers

Freeze/Unfreeze

The Freeze/Unfreeze functions are toggled via the console button labeled Freeze and the imaging screen button, labeled Freeze or Unfreeze depending on the mode.

When in Freeze mode there is no ultrasound imaging stream, and the Freeze/Unfreeze imaging screen button is labeled Unfreeze. Scrolling mode is automatically activated. Press the left trackball button to toggle scrolling mode on and off.

The Freeze function halts imaging on the last frame viewed and internally saves the last 50 frames imaged. The 50-frame buffer enables you to scroll backward to review previous frames and select the best one captured around the time the freeze was requested. When scrolling, the fps rate displayed onscreen indicates the speed at which you are reviewing frames.

Use the trackball to scroll through the buffered frames and the playback buttons to click through the sequence one frame at a time.

When in Unfreeze mode, there is a continuous ultrasound imaging stream in real time. The Freeze/Unfreeze imaging screen button is labeled Freeze.

Steering

The steering feature enables you to pan 26 degrees to the left or right without physically moving the probe. Steering can be helpful when:

- 1. A structure of interest is only partially visible, but cannot be viewed in full by physically moving the probe.
- 2. A structure is too large to be viewed in its entirety within the Zura system's 70 degree sector.

Steering is available from the keyboard during real-time imaging.

- Press ♡ > to pan right on image
- Press ♥ < to pan left on image

Steering also available from the trackball unless measuring or tracing mode is ON.

- Left-click on the right of the image to pan right
- Left-click on the left of the image to pan left

The steering angle advances in six steps to ± 26 degrees.

When steering is activated, two visual aids appear on screen (see Figure 50):

- 1. A set of blue arrows pointing in the direction of motion.
- 2. Blue tick marks beneath the image. One mark is positioned on the centerline of the image and is stationary; the other indicates the incremental movement of the steering.

To reset steering to the center position—i.e., angle = 0—press \bigcirc /.



Fig. 50 Single View real-time screen with steering visible

Image Acquisition

Press the Acquire button to save real-time images and ECG data as a three-second cineloop.

The B-mode acquisition frame rate accelerates from 24 fps to 50 fps.

The CF acquisition frame rate is twice the CF frame rate in continuous imaging up to a maximum of 30 fps. For example, if the CF frame rate is 13 fps, the acquisition rate is 26 fps. However, if the CF frame rate is 20 fps, acquisition is 30 fps. Figure below shows the Single View imaging screen during acquisition.

Image acquisition can be accomplished in one of two ways:

- 1. pressing the Acquire key on the console
- 2. pressing the Acquire button on the probe handle

After the three-second cineloop is captured it is written to disk. Figure 51 below shows the progress bar displayed when the cineloop is saved.



Fig. 51 Cineloop save progress bar

While the cineloop is being saved to disk, you may enter a comment in the cineloop comment window. The Comment window appears as soon as the Zura system begins saving the cineloop (see Figure 52). Comments are visible onscreen when previewing cineloops and when a cineloop is loaded.

You can set the hTEE view type by clicking one of the four view radio buttons. The view type is displayed in the cineloop load dialog and allows for quick cineloop identification. It is also used by the system to select automatically the relevant measurement package.

Entering comments is optional; comments can be added or edited later. If you don't wish to enter a comment, click the OK or Cancel button to clear the dialog box. Pressing the Acquire button on the probe handle or the console will also clear the dialog box.



Fig. 52 Cineloop comment dialog box

By default, the new cincloop is automatically loaded and ready for playback after acquisition. This behavior can be changed in the Acquisitions tab of the Configuration window. For additional information, see "Acquisition," page 80.

The Zura system stores up to 100 cineloops per patient. If a cineloop is acquired after the limit has been reached, you will be prompted to overwrite the oldest cineloop with the new acquisition. If desired, you can cancel the operation and manually delete selected cineloops. For additional information, see "Deleting a Cineloop," page 62.

Safety Feature

Recovery mode allows the ClariTEE probe to cool.

Recovery Mode

While the cineloop is being written to disk, the ImaCor Zura system software enters a nine-second recovery mode to allow the probe to cool. During the recovery period the acquisition function is disabled. The recovery frame rate is 6 fps.

Imaging is allowed after the cineloop has been saved to disk. At the end of the recovery period, acquisition is re-enabled and imaging automatically returns to the default frame rate. For B-mode imaging, the default frame rate is 24 fps. For CF imaging, the frame rate is inversely proportional to the size of the color box; i.e., the frame rate decreases as the size of the color box increases.

Loading a Cineloop

- 1 If you are starting from the Patient Information screen:
 - Select a patient name from the patient list
 - Click the Next button to enter the imaging environment
- **2** From the imaging environment, click the Load button. The Load Loop window is displayed.



Fig. 53 Single View mode – load cineloop dialog

- **3** Folders containing cineloops are listed in the Select Loop section. The folders are labeled with the date the cineloops were acquired. Locate the folder containing the cineloop you wish to view.
- **4** Double click the appropriate folder to see a list of cineloops contained within it. Each cineloop name is preceded by an icon that indicates whether it is a B-mode or Color Flow cineloop, as follows.
 - B mode icon:
 - CF mode icon:

To assist with cineloop identification, the first 32 characters of the cineloop comment are displayed.

Cineloops are also labelled with the date and time of acquisition. Locate the cineloop you wish to view.

- 5 Click once on the desired cineloop. If the Show Preview option is checked, the cineloop is displayed in the lower right of the Load Loop dialog window, along with comments and LVEDA, LVESA and FAC data, if any. See Figure 53.
- 6 Click the OK button to load the selected cineloop.

Note: Cineloops are automatically loaded after acquisition. This behavior can be changed in the Acquisitions tab of the Configuration window. For additional information, see "Acquisition," page 80.

Reviewing a Cineloop

Cineloop Image Quality Control

Cineloop image quality controls enhance existing loops for review purposes. Changes to the image control settings are not saved with the cineloop.

BRIGHTNESS CONTROL

The brightness control knob is located on the console. It allows you to increase or decrease the brightness of your image.



Fig. 54 Brightness control display

DEPTH BRIGHTNESS COMPENSATION (DBC)

DBC imitates TGC and applies brightness curves to the image.



Fig. 55 Depth brightness control (DBC) display

CONTRAST CONTROL

Contrast increases or decreases the contrast between light and dark areas in an image.



Fig. 56 Contrast control display

DEPTH CONTROL

Depth controls how much of an image you wish visible.



Fig. 57 Depth control display

Cineloop Control Buttons

Once the cineloop is loaded, a set of cineloop playback controls appears on screen, as shown in Figure 55.

Table 8 Cineloop Control Buttons

Button icon	Button name	Function
	Play/Pause toggle	Plays a cineloop.
1	Play/Pause toggle	Pauses a cineloop.
	Forward	 While a cineloop is paused, the Forward button advances the cineloop one frame.
		 While a cineloop is playing, the Forward button increases the playback speed. Playback speeds are:
		0.25 0.5 0.75 1.0 1.5 2.0
	Backward	 While a cineloop is paused, the Backward button reverses the cineloop one frame.
		 While a cineloop is playing, the Backward button decreases the playback speed. Playback speeds are: 0.25 0.5 0.75 1.0 1.5 2.0
(Jump	Jump to the next ECG marker

Figure 58 below shows where the cineloop controls are located on the imaging screen.



Fig. 58 Single View mode with cineloop control buttons

SCROLLING MODE

Scrolling mode enables you to scroll through frames instead of clicking frame by frame. Scrolling mode acts as a toggle; left click on the image to enable and disable scrolling.

- Roll the trackball right to scroll forward through frames
- Roll the trackball left to scroll backward through frames

COLOR FLOW INDICATOR ICON

When Color Flow (CF) images are available for a cineloop, the CF indicator icon is displayed. Click the CF indicator icon to toggle on and off the display of Color Flow images during cineloop playback and when viewing frozen real-time images.

Table 9 Color Flow Indicator Icons

lcon	Function
CF	Indicates that color flow images are available and are displayed.
CF	Indicates that Color Flow images are available but are not displayed.

Exporting a Cineloop as a Movie

Cineloops can be exported as movies that can be viewed on standard video players. Exported movies must be saved to a USB thumb drive.

To export a cineloop as a movie:

- 1. Press the Load button to load the cineloop.
- 2. Insert a USB thumb drive.
- Right-click to view the global context menu and select the Export Movie option.
- 4. When the Export Movie window opens, specify the following;
 - Desired movie format (AVI or MOV)
 - Scaling, if desired
 - Inclusion of ECG trace, if desired
- 5. Click the Export button (see Figure 59). A progress bar is displayed while the movie is exporting (see Figure 60).

Saving the Current Frame as a BMP Image

The current frame can also be saved as a BMP image, as shown in Figure 59. Follow the steps for saving a cineloop as a movie, but select the BMP (single-frame) option from the Save As window. Like cineloops, single-frame BMP images are also watermarked with the ImaCor logo.



Fig. 59 Export movie format options

Export Movie	×
Select the movie format and AVI MDV	d options: Scale[0.1-1.0] 1 Include ECG 🔽 Whole Screen 🗖
Exporting the movie Pleas	se wait.
	Cancel

Fig. 60 Export movie progress bar

Measurement Tools

The Zura system software provides two measuring tools:

- 1. The Tracing Tool, for area measurements. This tool is key when performing hTEE calculations.
 - FAC
 - RV/LV ratio
 - SVC collapsibility index

2. The Measuring Tape Tool, for distance measurements.

The Tracing Tool

To access Tracing Tool, press the Trace key on the console or click the onscreen Measure button.

When using the FAC Wizard to guide an FAC calculation, the Tracing tool is automatically activated. When trace mode is ON, the backlit key color changes from the default blue to green. In addition, the informational message "Trace: ON" is displayed on the screen.



Fig. 61 Accessing the Tracing tool with the on-screen Measure button.

CREATING A TRACE

Right click the trackball to set the first anchor point and begin the trace. Once you start tracing, any movement of the trackball sets points at a fixed interval. The interval is called trace resolution and can be changed in the Configuration dialog box (see page 79).

Traces are automatically saved with the patient exam.

A new trace can be used to update a previous trace without using the FAC wizard. To substitute a new LVEDA or LVESA value:



- 1. Right click the desired field while Tracing or Measuring mode is on and a closed trace is displayed.
- 2. Select the "Use current trace" option.

LVEDA:	14.6 cm ²
LVESA:	7.1 Use current trace
FAC:	51 %
Date:	Nov-12-2007
Time:	01:01:17 PM

Fig. 62 Manually updating a trace value

The date and time the cineloop was acquired is also displayed in the FAC information area.



Fig. 63 Area measurement – 12 cm depth

When tracing, the trackball is set to a lower sensitivity to reduce tracing errors. Trackball speed during tracing can be customized in the Configuration dialog box (see page 79). A low speed number indicates low sensitivity; high speed number high sensitivity.

To backtrack one point while tracing, press the Delete or Backspace key. Pressing Delete multiple times erases multiple points.



Fig. 64 Trace with multiple points deleted

Right click to close a trace.

Turning off tracing mode hides the trace until tracing mode is toggled back on; it does not delete the trace.

You can have only one trace onscreen at a time. If you wish to create a new measurement, you must delete the existing trace.

EDITING A TRACE

When the tracing cursor is moved about an anchor point, the anchor point is highlighted and can be edited.

- Press the Delete key to erase the highlighted anchor point. This opens the trace at the location of the deleted point.
- Right click on an anchor point to open the trace at that location without deleting the point.
- Highlight multiple anchor points by holding down the Shift key and moving the trackball left or right; you need not follow the contour of the trace. Press the Delete key to erase the selected segment of the trace. You can then redraw the trace segment.

CHAPTER 4: IMACOR ULTRASOUND IMAGING SOFTWARE



Fig. 65 Single View mode, editing multiple points

DELETING A TRACE

You must be in Trace mode to delete a trace. Be sure no points are selected, then press the Delete key.

To clear all distance and area measurements at once, press the \bigcirc key and the Delete key in combination.

The Measuring Tape Tool

The measuring tape tool can be accessed in two ways: press the Measure console key or click the onscreen Measure button and select Distance. When measure mode is ON, the backlit key color changes from the default blue to green. In addition, the informational message "Measure: ON" is displayed on the screen.





Fig. 66 Accessing the Measuring Tape tool with the on-screen Measure button.

TAKING A MEASUREMENT

Right click the trackball to set the first end point and begin measuring. Right click again to set the second end point, completing the measurement.

To create a second measurement, right click to set a new starting point, and right click again to set the second end point. You can set a maximum of two distance measurements, which will be labeled on screen as Dist1 and Dist2. The end points are labeled 1 and 2 (see Figure 68).

If the end point of a measurement must begin on top of or very near another point, move the cursor to the desired end point location and pause until the highlight on the neighboring point disappears. Then right click to set the point in the new position.

Distance measurements are automatically saved with the patient exam.



Fig. 67 Single View mode, taking area measurements (Trace: ON)

EDITING A MEASUREMENT

Move the tape measure tool about one of the end points (caliper). When the end point is highlighted, right click. This frees the end point, which becomes a floating point. Move the floating point to the desired position, then right click again to reset the end point.

DELETING A MEASUREMENT

Move the tape measure tool about either of the end points. When an end point is highlighted, press the Delete key.

To clear all distance and area measurements at once, press the \bigcirc key and the Delete key in combination.

Viewing Measurements

You can view distance and area measurements at the same time. Select both modes, and both measurements are displayed simultaneously.



Fig. 68 Single View mode with distance measurements (Measure: ON)

If you turn off either trace or measure mode, both measurements are hidden. Reactivating one mode will display only the measurement associated with it. The second measurement can be turned on separately, if desired.

VIEWING MEASUREMENTS ON CINELOOPS

When a cineloop containing distance and/or area measurements is loaded, the measurements are visible only if measure mode is ON. *The measurements displayed may not apply to the current frame.* To view the frame on which the measurements were based:

- Move your mouse over the distance or area value. A blue rectangle appears around the value.
- Left click the measurement value to jump to the corresponding frame.

Fractional Area Change (FAC) Calculations

Fractional area change (FAC) calculations are based on measurements of the left ventricular end-diastolic area (LVEDA) and the left ventricular end-systolic area (LVESA). The ImaCor Zura system software provides two ways to complete an FAC calculation:

- 1. The FAC Wizard, which guides you through the process of completing an FAC calculation. The FAC Wizard is available only in Single View mode.
- 2. Direct measurement by saving a current trace as the LVEDA or LVESA (see "Creating a Trace," page 51).

The FAC Wizard

Using the FAC Wizard to complete a calculation is a six-step process.

1 Click the FAC Calc button onscreen to begin. If you have not loaded a cineloop, a message instructing you to do so is displayed.



Fig. 69 Load cineloop instructional message

2 With a cineloop loaded, step or scroll through the loop and select the end-diastolic frame.



Fig. 70 FAC calculation in Single View mode – end diastolic identification and trace

3 Right click the trackball to set the first anchor point, then begin tracing. When you complete the trace, right click again to close the trace. The LVEDA value is displayed on screen, as shown in Figure 71.



Fig. 71 FAC calculation in Single View mode – end systolic identification and trace

- **4** Click the ES button. Scroll or step through the frames and select an end-systolic frame.
- **5** Trace the LVESA.
- **6** When you have finished, click the Done button. The LVEDA, LVESA, and FAC values are automatically calculated and displayed, as shown in Figure 72.



Fig. 72 Completed FAC calculation

FAC data are saved with the loop and are displayed when the loop is previewed and loaded. Cineloop acquisition date and time are displayed under FAC measurements.

REVIEWING A TRACE

To review an existing FAC trace, click on the desired measurement (LVEDA or LVESA). The Zura system software automatically displays the traced frame and highlights the trace value.



Fig. 73 FAC calculation results

TRACE OVERLAY MODE

Trace overlay mode allows you to display or hide the LVEDA and LVESA traces on cineloops. When overlay mode is active, the color of the backlit Overlay console key changes from the default blue to green.

Trace overlay mode can be turned **on** three ways:

- 1. Activate the FAC Wizard to automatically turn on trace overlay mode.
- 2. Left click the LVEDA or LVESA value.
- 3. Press the Overlay console key to toggle overlay mode on and off.

Trace overlay mode can be turned **off** by pressing the Overlay console key.

SUPERIMPOSED COMPARISON

When two cineloops are loaded in Split View mode and both have FAC calculations, traces can be compared by superimposing one on the other. See page 62 for details.

Deleting a Cineloop

You can store a maximum of 100 cineloops per patient. As you near a patient's cineloop limit, warning messages are displayed, as shown in Figure 74.



Fig. 74 Warning: cineloop limit approaching

Once 100 cineloops have been saved, you cannot acquire additional images unless you overwrite the oldest loop or delete one or more cineloops. For additional information, see the master message listing, page 101.



Fig. 75 Warning: cineloop limit reached

To delete a cineloop:

- 1 Press the Load button.
- **2** From within the Load Loop dialog screen, select a cineloop and right click. Multiple cineloops may be selected for deletion using the Shift or CTRL key.
- **3** Press the DEL or Backspace key, or select the Delete option when the popup menu appears.



Fig. 76 Delete cineloop screen

4 When the warning message appears, click the Yes button to delete the cineloop and No to cancel the operation. Once a cineloop has been deleted, it cannot be retrieved.

Warning	1	×
	Patient: (A5C40655-9084-4ECE-A12A-7A46F34DBADA) Do you want to permanently delete the selected patient data?	
	Yes No	

Fig. 77 Delete cineloop confirmation

Split View Mode

Split View mode enables you to evaluate side by side either two cineloops (Figure 78) or a cineloop and real-time imaging (Figure 80). The Views key on your console acts as a toggle between Single and Split View.

Real-time imaging is always displayed on the right of the screen. The active loop name is always highlighted.

Steering is available in Split View mode provided real-time imaging is taking place.



Fig. 78 Split View mode – cineloop comparison

Loading Two Cineloops in Split View Mode

- 1 Click the Load button.
- **2** When the Load Loop window is displayed, select a loop from the left and the right pane, as shown in Figure 79.
- **3** When the cineloops have been selected, press the OK button.

elect Left Loop:		Select Right Loop:	
04-15-2011 04-14-2011		04-15-2011 04-14-2011	
20110414-224952		20110415-080214	
20110414-192841	SAM	20110415-080137	
20110414-180000		20110415-080030	
20110414-175825		20110415-075958 SAM	
20110414-152441	TGSAV	20110415-055832	
20110414-152354		20110415-025306	
20110414-145027	SAM	20110415-025140 TGSAV	
20110414-144825		20110415-024759	
Show Preview	and the second sec	☑ Show Preview	
GSAV		TGSAV	
election: 0110414-152441		Selection: 20110415-025140	

Fig. 79 Split View mode - loading two cineloops

Viewing a Cineloop and Real-Time Images in Split View Mode

- 1 Click the Load button.
- 2 When the Load Loop window is displayed, select a cineloop from the left pane.
- **3** When the cineloop has been selected, press the OK button. The cineloop is displayed on the left of the Split View screen, real-time imaging on the right.

Note: If you load a cineloop from the right pane during real-time imaging, the cineloop replaces the real-time view. Click the Unfreeze button to resume real-time imaging.

CHAPTER 4: IMACOR ULTRASOUND IMAGING SOFTWARE



Fig. 80 Split View mode – real-time image and cineloop comparison

Synchronized Playback

The synchronized playback function enables you to review a full cardiac cycle from two different cineloops in a synchronized fashion. You must view cineloops—the function is unavailable during real-time imaging—and each loop must have at least two ECG markers.

The software synchronizes the playback using the first two consecutive ECG markers it encounters on each cineloop. If the cineloop frames were not advanced after being loaded, these will be the first two ECG markers in the cineloop.

If one or both cineloops contain more than two ECG markers, you can specify which marker should be used to begin the synchronization.

- **1** To advance to the next ECG marker, click on the Jump to ECG marker button, located below the corresponding cineloop.
- 2 Repeat the process on the other cineloop, if desired.



Fig. 81 Synchronized playback

ACTIVATING AND DEACTIVATING SYNCHRONIZED PLAYBACK

The Sync button, shown in Figure 82, toggles synchronized playback on and off.



Fig. 82 Playback buttons with Sync button active.

PLAYBACK SPEED CONTROLS

During playback, the Forward and Back buttons act as speed controls. Forward increases the frame rate and Back slows the playback. Playback speeds are:

0.25 0.5 0.75 1.0 1.5 2.0.

To view a frame-by-frame synchronized playback, you must scroll on the master (right) side of the display.

When in pause mode, the Forward and Back buttons enable you to step forward and backward one frame at a time.

Superimposed Mode

Traces can be superimposed for comparison. The Compare ED/ES menu is enabled when two cineloops are loaded and each cineloop has FAC data.

To access the menu, right click on the image and select Compare ED/ES from the popup, as shown in Figure 83.



Fig. 83 Compare ED/ES menu option

The Superimposed View window is displayed in the center of the screen. By default, the LVEDA trace of each loop is superimposed on a black background. The LVEDA trace from the left panel is in green, and the right in blue (see Figure 84).



Fig. 84 Split View mode with Compare ED/ES selected

To view superimposed LVESA traces, click the ES button on the Superimposed View window.

Color Flow Imaging

Color flow (CF) imaging, also known as Color Doppler, uses color to detect blood flow and determine flow direction. Colors in the middle to the top of the Color Map indicate positive blood-flow velocities; i.e., blood flow toward the transducer. Colors in the middle to the bottom of the range show negative blood-flow velocities; i.e., blood flow away from the transducer.

Color Flow mode can be activated in real-time imaging mode. You can activate CF imaging mode by pressing the "C" key on the console or by clicking the Color Flow button in the button bar.

When Color Flow imaging is on the "Color ON" indicator is displayed, as shown in Figure 85.



Fig. 85 Real-time imaging with Color Flow imaging active

The Color Box

The color box, shown in Figure 86, defines the region where CF data is displayed. The color box is visible when viewing:

- real-time images
- frozen buffer images
- cine playback images

The color box is not displayed when CF overlay is turned off.

When Color Flow imaging is activated, the color box appears in edit mode (see Figure 86). There are two edit modes for the color box:

- 1. Edit Mode 1 allows you to position the center of a fixed-size color box by moving the trackball horizontally and vertically.
- 2. Edit Mode 2 allows you to resize the color box by moving the trackball horizontally (width) and vertically (height). The color box's center does not move when it is being resized (see Figure 87). Edit mode tips are displayed below the Real Time label.

Left click on the ultrasound sector to toggle between edit modes. After changing the color box position or size, left click for resultant CF overlay changes to immediately take effect. Alternatively, wait 1.5 seconds and the revised CF overlay will take effect automatically.

Edits to the size and position of the color box also take effect in 1.5 seconds

Color flow image quality controls are explained in "Real-Time Image Quality Controls for Color Flow Imaging," page 40 of this document.



Fig. 86 Color Flow imaging with color box in edit mode
CHAPTER 4: IMACOR ULTRASOUND IMAGING SOFTWARE



Fig. 87 Sizing the color box

B-Mode Priority

B-mode priority is a user-adjustable grayscale value between 10 and 255. The Bmode data takes priority when the grayscale value is greater (brighter) than the user-set priority value. For example, if you set a priority level of 150, B-mode data is shown instead of CF data when the B-mode data value is greater than 150. This prevents the display of CF data with no clinical benefit, such as in tissue.

Notes:

- 1. B-mode data values range from 0 (black) to 255 (white).
- 2. Priority is also available on cineloop playback.

Like other image quality parameters, the priority setting is saved as a part of the patient record and restored when the patient is again imaged.

When a cineloop is acquired, it is stored with no B-Mode priority. You can adjust the priority on playback.

In split view, Priority can only be adjusted on the right pane cineloop.

SETTING B-MODE PRIORITY

To set B-mode priority, click the grayscale gradient bar on the right of the screen. A red horizontal bar appears indicating the position on the scale where priority has been set, as shown in Figure 88. Point and click on the scale to change the location of the bar and priority setting.



Fig. 88 Grayscale gradient bar showing priority setting

Velocities

The velocities menu is displayed by right or left clicking the velocity display (V Max) during Color Flow imaging.

There are three velocity modes to choose from: Low, Medium, and High.

CHAPTER 4: IMACOR ULTRASOUND IMAGING SOFTWARE



Fig. 89 Velocities menu with color box between 9 and 11 cm

Velocity setting availability is associated with the depth of the color box. When the bottom edge of the color box exceeds 9 cm, the velocity setting automatically changes from High to Medium. When the image depth exceeds 11 cm, the velocity settings automatically shifts from Medium to Low.

To manually select a velocity setting, right-click the velocity display and choose from the velocity pop-up menu.

Note: In low velocity mode, the velocity text display flashes as a reminder that you are operating in low velocity mode, which is not the primary velocity mode.

Color Maps

Five Color Maps are available. To switch from one Color Map to the next simply left-click on the Color Map.



Fig. 90 Highlighted color map

CF Overlay in Playback Mode

When in frozen buffer mode, or in playback mode, Color Flow images can be shown or hidden by clicking on the CF icon above the images.



Fig. 91 CF images displayed, Color Flow



Fig. 92 CF playback in split screen mode.

Cineloops acquired with Color Flow mode ON have a distinct icon so that they can be distinguished from B-mode cineloops in the Cineloop Load Dialogs.

elect Left Loop:	Select Right Loop:
C 03-11-2010 C 02-08-2010	03-11-2010 02-08-2010
20100311-123617 TGSAV	20100311-123617 TGSAV
20100311-123544 4 chamber view	20100311-123544 4 chamber view : 30
20100311-104409 cf gain: 24 cf filt	20100311-104409 cf gain: 24 cf fill
20100311-104349 Crigan: 24 crint	20100311-104349 Cligan: 2 <u>4</u> Clint
cf gain: 24 cf filter: 14	cf gain: 24 cf filter: 14
Selection:	Selection:

Fig. 93 CF-mode and B-mode cineloops have distinct icons

Printing

A printout contains the following information:

- Patient first and last name
- Patient ID
- B-mode image
- CF image (in grayscale because the printer is black-and-white only)
- FAC data, if present
- Distance measurements, if present
- Area measurement, if present
- System date and time
- Cineloop date and time, if printing a cineloop
- Cineloop comment, if present
- Institution name

Single View Mode

Press the Print button to produce a hardcopy of the active frame, as well as FAC measurements and comments, if applicable. The ECG waveform is not included in the printout.

Split View Mode

Press the Print button and select from the print dialog box the image you wish to print.



Fig. 94 Print dialog box in Split View mode



Fig. 95 Sample printed output

The warning message below is displayed if you attempt to use a printer that isn't properly connected to the system. For more information, see the master message listing on page 101.



Fig. 96 Printer connection error message

Configuring the System

The Configuration dialog window consists of seven tabbed sections:

- General
- Acquisition
- Date-Time
- DICOM
- Reports
- Language
- Version
- Help

The Configuration dialog can be accessed from either the Patient Information screen or the imaging environment. Imaging-related options can be configured only from the imaging environment.

General

The General tab enables you to change the Imaging and Tracing/Measuring default parameters.



Fig. 97 Configuration > General tab

IMAGING DEFAULTS

After changing Imaging defaults, restart the system to apply the new settings.

MEASURING AND TRACING SETTINGS

Trace Resolution

Trace resolution is the interval between two consecutive anchor points. This interval is expressed in pixels, so a smaller resolution value results in a larger number of anchor points. The range is 1–20.

Trace Trackball Speed

When tracing is enabled, the speed of the trackball is automatically lessened to reduce errors. The Trace Trackball Speed option enables you to specify a custom trace speed. The speed range is 1–20, with 1 being the lowest speed and 20 the highest.

Trace Color

The Trace Color option enables you to specify anchor point color.

Measure Color

The Measure Color option enables you to specify the color of measurement lines drawn with the Tape Measure tool.

INSTITUTION

Enter the name of your institution (hospital, medical center, etc.) in this field. Your institution name will be displayed on all screens. Limited to 30 characters.

RESET BUTTON

Press the Reset button to restore the manufacturer's default values. If you press the Reset button, the confirmation shown below is displayed.



Fig. 98 Reset button confirmation screen

Acquisition

The acquisition tab lets you select one of the following cineloop acquisition options.

- The "Load acquired cineloop" option directs the system to load a newly acquired cineloop after it has been acquired. This is the default option.
- The "Unfreeze to real-time" option returns the system to real-time operation after a cineloop has been acquired.

Selecting the "Popup comment dialog" option directs the system to display a cineloop comment dialog pop-up after cineloop acquisition. This is the default setting; it can be specified with either of the above options.



Fig. 99 Configuration >Acquisition tab

Date-Time

The Date-Time tab enables you to set the current date and time.

Changing the system date and time has no effect on probe life.



Fig. 100 Configuration > Date-Time tab

DICOM

Digital Imaging and Communications in Medicine (DICOM) is the standard format for distributing and viewing all types of medical images. DICOM files contain the following information:

- Patient information
- Cineloop
- ECG
- LVEDA, LVESA, and FAC, if present

The DICOM tab in the Configuration menu enables you to:

- Specify that cineloops be exported in DICOM format
- Indicate if JPEG compression should be applied
- Configure a DICOM Archiving Server (application entity title [AET], Host, Port)
- Configure a DICOM MWL Server (AET, Host, Port)
- Customize DICOM synchronization rules

- Customize system network properties
- Edit the client application entity title (AET)

Configuration	×
General Acquisition Date-Time DICOM L	anguage Reports Version Help
🌄 USB Export	🗾 Archiving Server
Export DICOM Images	Active AET
Cverwrite DICOMDIR	Host
Compression	Port Test
🔽 Use Lossy JPEG Compression	Rules PatientID
📝 MWL Server	PatientDOB
Active AET	PatientSex
Host	Client AET ZURA01
Port Test	Network
	OK Cancel

Fig. 101 Configuration > DICOM tab

Patient Export Options

If the Export DICOM option is selected, DICOM files are automatically created as cineloops are exported to the DTM. If Color Flow images are available, they are embedded in the DICOM file.

If the Overwrite DICOMDIR when Exporting option is checked, the DICOMDIR file is overwritten with each patient export.

Notes:

- If multiple patient records are selected for export and the Overwrite DICOM-DIR when Exporting option is checked, only the last record selected will appear on the DICOMDIR.
- DICOMDIR and DICOM files are encrypted before they are saved onto a DTM. An ImaCor decryption application is provided to import DTM data into a DICOM workstation.

Archiving Server Configuration

To configure the Archiving Server, the "Active" box must be checked. The server configuration parameters are:

- Target AET
- Host name / IP Address
- Port number

By default the client AET for Zura is ZURA. It can be changed by entering an appropriate name in the Client AET field .

The Network configuration parameters are those exposed by the Windows XP operating system and will require an ImaCor technician (locally or with remote access) to set up.

Once the server has been configured the user can click on the "Test" button to check communication with the server.

Checking the Archiving Server "Active" checkbox will enable the "Sync" button on the Patient Screen.

Custom Synchronization Rules

The Zura System provides six basic rules that can be turned ON or OFF independently to customize the synchronization process. These rules are:

- 1. Require Cineloop Comment
- 2. Require Patient ID
- 3. Require Patient DOB
- 4. Require Patient Accession Number
- 5. Require Patient Sex
- 6. Max Daily Cine

For example if rules (1) and (2) are turned ON, only cineloops that have a comment and are for a patient whose ID is not empty will be sent to the Archiving Server. The rules can only be turned ON/OFF by an ImaCor technician.

If the Require Patient Accession Number rule is ON, cineloops are sent to the server only if the accession number field in the patient information screen is completed for the corresponding patient.

The Max Daily Cine rule enables you to specify the maximum number of cineloops in a daily folder that may be sent to the archiving server.

Synchronization operations occur when the Sync button is clicked. A synchronization log provides information about the last Synchronization operation that was performed. To view the Log, simply right-click on the "Sync" button and choose "Show Synchronization Log."

Modality Worklist (MWL) Server Configuration

To configure the MWL Server, the "Active" box must be checked. The MWL Server configuration parameters are similar to those required by an Archiving Server.

Once the server has been configured the user can click on the "Test" button to check communication with the server.

Checking the MWL Server "Active" checkbox will enable the "Check Worklist" menu in the Patient list.

Compression

If the Compression option is selected, JPEG compression is automatically applied to the DICOM files. The compression option applies to both exported data and data sent to the DICOM Archiving Server (if configured).

Language

The Language tab enables you to select from a drop-down menu the language used in the software interface. Only languages appropriate for your location will be present on the menu.

onfigura	tion				
General	Acquisition	Date-Time DICOM	Language	Reports Version	Help
Select	one of the in	stalled languages to u	se in the inter	face:	
				_	
Ei	nglish US			•	

Fig. 102 Configuration > Language tab

Reports

The Reports tab enables you to enable/disable the hTEE reporting feature.



Fig. 103 Configuration > Reports tab

HTEE REPORTS

When the hTEE report feature is enabled, you can record information about each exam/imaging session in the exam log. The exam log has three tabs:

- pre-hTEE
- post-hTEE
- summary

The pre-hTEE tab, shown in Figure 104, is where information is recorded about the indications for imaging and current management. The name of the physician performing the exam must be entered.

Exam Dates:		Exa	ims:			
01-22-2013			Start Time	End Time	1	Export on Txt
		11	12:02:12 PM	(Current)		
hehTEE Post-hTEE Summary Indications:			Current Management:			
🔽 Hemodynamic Instability	Possible Pericardial Collection	1	₩ IV Fluids		Е ЕСМО	
Fluid Management	Management of Cardiovascular Support	1	Vasopressors		IT LVAD	
Possible LV Dysfunction	C Other	1	Inothopes		T RVAD	
Possible RV Dysfunction		1	🏳 Pulmonary Vasodil	stors	□ Other	
Notes:		Ð	Physician			
[1	Dr Who		1	

Fig. 104 Exam log – pre-hTEE tab

The post-hTEE tab, shown in Figure 106, is where information is recorded about views obtained, hTEE assessment, and hemodynamic management.

Exam Dates;		Exams:				
01-22-2013		Start Time	End Time (Current)			Export os Txl
ehTEE PosthTEE Summay	UTEE Assessment	1		Heredree	nic Managana	-t-
views obtained.	TILE ASSESSMENT.			No Chapon	Increased	Decreased
TIGSAV E ME4CH	F Hypovolemic	esponsive I Collection	IV Fluids	Г	प्र प	Г
	RV Dysfunction		Vasopressors	Г	Г	9
	Comment:		Inotropes	Г	Г	Г
Notes:	F	_	Pulm Vasodilators	Г	Г	Е
			ECMO, LVAD, RVAD	Г		Г
			-	- c	E.	E

Fig. 105 Exam log – post-hTEE tab

		E xams.
01-22-2013		Start Time Cod Time Export As T at I1:05:06 AM (Current)
Yes	No	
П	Г	Describe how hemodynamic management was individualized\nand guided over time:
Г	Г	
г	Г	
г	Г	
	Yes - -	Yes No r r c c c c

After the final exam is done, the summary report must be completed.

Fig. 106 Exam log – Summary tab

Missing information in any tab is highlighted in red.

Exam reports are stored by date in daily folders, which are visible in the the Exam Date field. When a folder is selected, exams completed on that date are in the Exams field, along with the starting and ending times of the exam.

Click the Export as Txt icon in the upper right of the Reports dialog box to export hTEE report data to a USB drive. The report is titled with the patient's GUID followed by "-Report.txt."

Version

The Version tab enables you to check the version of the Zura software installed on the system.



Fig. 107 Configuration > Version tab

Software Update

Upgrade packages are machine-specific and require a registration code that is displayed in the Configuration Dialog >Version tab (see Figure 107).

DOWNLOAD THE SOFTWARE UPGRADE AND REQUEST A LICENSE FILE

- 1 Log into the ImaCor support web page: http://www.imacorinc.com/support.html
- 2 Click on the software upgrade you wish to download.
- **3** Enter the Zura machine registration code, which can found on the Version tab of the Zura Configuration window.
- 4 Click on the Go button. A license request is automatically sent to the ImaCor support team.
- **5** Click on the "Download File" option.
- 6 Navigate to the location on your computer where you wish the file saved.

UPGRADE THE ZURA SOFTWARE

- 1 Save the .icm file that was emailed to you onto a USB thumb drive. For example, if your thumb drive is the F: drive, save the .icm file to F:
- **2** Unzip the software downloaded upgrade and move the .exe file inside the zip folder to your USB thumbdrive. The .icm and .exe files must be placed at the same location.
- 3 Plug your USB key into the Zura system you are upgrading.
- **4** Go to Configuration > Version tab and click the Update button. You can accomplish this from either the Patient Information screen or the Imaging screen.
- 5 Select the .icm file on your thumbdrive and click "OK."
- **6** When the Upgrade Wizard appears (it may take up to 5 seconds), follow the instructions.
- **7** Do not disconnect the USB thumb drive until a message is displayed instructing you to reboot the system.
- **8** After the system reboots you will be prompted by the Upgrade Wizard to click on "Finish." Once you've clicked Finish, the upgrade is complete.

Help

In addition to the online Help system (see page 19), the Zura system provides three other user assistance options:

- 1. ToolTips
- 2. Remote Assistance
- 3. Analysis Package



Fig. 108 Configuration > Help tab

Tool Tips

If the Show ToolTips option is selected, informational pop-ups are displayed when you hover the mouse over interactive features like software buttons and dialog box elements.



Fig. 109 ToolTips visible in imaging screen

Remote Assistance

The Remote Assistance button enables you to establish a link with an ImaCor technician who can access the Zura system remotely for troubleshooting.

Once the Remote Assistance button is clicked, the Support Internet Explorer window is displayed.

- Enter the 6-digit PIN code provided by ImaCor tech support.
- Click the "Connect to technician" option to establish a secure connection with a technician's console.

🖉 Support - Windows Internet Explorer	×
ImaCor	×
Enter your 6-digit PIN code:	
Connect to technician	
	Y

Configuration > Support tab

Analysis Package

The support technician may request an analysis package, which is a file that provides information about your system and assists in troubleshooting.

To create an analysis package:

- Select Configuration > Help
- Insert an external USB drive into the USB port in the Zura system. You need not use the special ImaCor data transfer module; any operational USB drive will work.
- Click the Analysis Package button

Email the analysis package (ZuraPkg.zip) to support@imacorinc.com.

If there is insufficient free space on the USB drive to accommodate the analysis package, or if the USB drive is removed prematurely, the error message shown in Figure 110 is displayed.

- In the case of insufficient free space on the USB drive, you may:
 - Remove the drive from the Zura system.
 - Insert the drive in any personal computer with a USB port.
 - Delete unnecessary files. Files may not be deleted from the Zura system.
 - Retry the operation
- If the USB drive was removed prematurely i.e., before the analysis package is completely written to the drive you may:
 - Reinsert the drive in the Zura system.
 - Retry the operation.

The error code shown in the message may vary and is for the use of ImaCor technicians only.



Fig. 110 Analysis package export error

If no USB drive is present when the Analysis Package button is pressed, the error message below is displayed.



Chapter 5: Safety

Read this safety information before using the ImaCor Zura systems. The information applies to the entire system: machine, umbilical, disposable probe, and accessories.

Disposable Single-Use ClariTEE

The ClariTEE is a single-use probe provided in sterile packs.

During episodic assessments, the probe tip is flexed upward via the articulation lever on the handle. Probe pressure on mucosal tissue during flexing requires continuous operator action; there is no locking mechanism in handle.

The indwelling portion of the ClariTEE probe is detachable from the control handle. Detaching the control handle allows for improved patient mobility and comfort during probe dwell time.

The distal end of the ClariTEE probe may be immersed in water or other gastric liquids. The ClariTEE probe's distal end and insertion tube meets the requirements of WATERTIGHT, (IPX7 rating) per ingress protection standard IEC 60529

(I) Caution

- The ClariTEE probe is designed for single patient use only. If the probe is reused there are known characteristics and technical factors that could pose safety risks as follows:
 - The probe is not designed to survive Cydex disinfection. Cydex disinfection of the ImaCor ClariTEE probe could present a toxicological and electric shock risk to the patient
 - The electrical insulating properties of the ClariTEE probe are unknown if the probe is reused. Reuse of the ClariTEE single-use probe could present an electric shock risk to the patient.
 - The biocompatibility properties of the ClariTEE probe are unknown if the probe is reused. Reuse of the ClariTEE single-use probe could present safety risks due to material incompatibility with the patient.
 - The mechanical integrity of the ClariTEE probe is unknown if the probe is reused. Reuse of the ClariTEE single-use probe could present a mechanical risk to the patient; i.e sharp edges.
- Flex the probe only when imaging.
- Disconnect the probe from the handle when not conducting an episodic assessment.
- The maximum dwell time for the ClariTEE probe is 72 hours.
- The ClariTEE probe is for single use only and is to be disposed after single use or if other wised contaminated by blood, body fluids, or other biological materials. Do not sterilize or disinfect the ClariTEE disposable probe.
- Handle the probe carefully.
- Do **not** drop the disposable probe.
- Do not pinch or kink the disposable probe.
- Do not use probe if cracked, broken, or damaged.
- Do **not** force the probe during intubation. If resistance is met during intubation, gently withdraw the probe.

1 Caution

The ImaCor Zura system should not be used beyond a maximum cumulative imaging time per patient of 6 hours.

A Caution

Disconnect probe from handle when not conducting an episodic assessment.

Important: To restore imaging, the user should toggle the Freeze/Unfreeze command button.

Reboot the computer if the system ceases to operate or appears frozen; i.e., does not respond to keystroke commands. • Separate disposed probe(s) from normal waste. Refer to local regulations.

The ALARA Principle

The Zura system has been designed in accordance with the principle of ALARA. ALARA is an acronym for an important principle in ultrasound energy protection and stands for "As Low As Reasonably Achievable." The aim is to minimize the amount of acoustic energy exposure to that which is necessary to accomplish the intended diagnostic use of the device.

Applying ALARA

The Zura system is a track 1 ultrasound imaging device. This means that acoustic output levels do not exceed FDA limits established in 1976. Further, the output power has been set as low as reasonable to produce an image in the area of interest. Every appropriate design effort has been made to use the minimum acoustic output required to obtain clinically adequate image.

ALARA Operator Control

CUMULATIVE IMAGING TIME

The Zura system is intended to be used episodically. It is not a continuous monitoring device. The device should be in imaging mode only while conducting episodic assessments of patient cardiac function.

Normal use should consist of a maximum of six episodic assessments over a 24hour period with maximum intubation time (the probe dwells within the patient) not to exceed 72 hours. The duration of a typical assessment will be in the range of 5–15 minutes. Hence, the typical cumulative imaging time for each patient is in the range of 1.5 hours to 4.5 hours during a 72-hour indwelling time period.

When not conducting an episodic assessment the device should not be imaging; no acoustic energy shall be delivered to the patient. In order to terminate imaging mode:

- Click Freeze button, or
- Exit patient application

MECHANICAL INDEX

The Mechanical Index (MI) is an index derived from the acoustic output of ultrasound diagnostic systems. The lower MI is achievable at the 7 MHz frequency setting. The operator is advised to use the 7 MHz setting where reasonable to achieve the desired image.

ALARA Automatic Controls

TWENTY MINUTE SOFTWARE INTERLOCK

If the operator leaves the Zura system in imaging mode while not conducting a bedside episodic assessment, the acoustic output automatically times out after 20 minutes. This limits the acoustic energy unnecessarily delivered to the patient.

WATCHDOG MONITOR

In the event of a software or operating system failure, the watchdog monitor limits burst mode duration to three seconds.

Transesophageal Echo Imaging (TEE) Potential Complications

1 Caution

TEE procedures shall be conducted by a qualified physician.

TEE is considered a safe procedure in that the chances of adverse events are low. In rare instances the following patient risks are associated with transesophageal echo:

- Dental injury
- Bleeding or tearing of the esophagus or stomach
- Difficulty swallowing after the procedure
- Respiratory distress

Electrical Safety

The Zura system complies with the following safety standards.

IEC 60601-1 (1988) Medical Electrical Equipment – Part 1: General Requirements for Safety.

IEC 60601-2-37 (2001) Medical Electrical Equipment – Particular Requirements for the Safety of Ultrasonic Medical Diagnostic Equipment.

- To reduce the risk of electric shock, do not open system enclosures.
- To reduce the risk of injury, do not operate the system in the presence of flammable gasses or anesthetics. Explosion can result.
- To reduce the risk of electrical shock, use only properly grounded equipment. Shock hazards exist if the power supply is not properly grounded. Grounding reliability can only be achieved when equipment is connected to a receptacle marked "Hospital Only" or "Hospital Grade" or the equivalent. The grounding wire must not be removed or defeated.
- To reduce the risk of electrical shock, inspect the transducer face, housing, and cable before using the transducer. Do not use the transducer if the transducer or cable is damaged.
- To reduce the risk of electrical shock, always disconnect the power supply from the system before cleaning the system.
- To reduce the risk of electrical shock and fire hazard, regularly inspect the power supply, AC power cord, and plug for damage.
- To reduce the risk of electrical shock, use only accessories and peripherals recommended or supplied with the system by ImaCor. Use of accessories and peripherals not recommended by ImaCor could result in electrical shock. Contact ImaCor or your local representative for a list of accessories and peripherals available from or recommend by ImaCor.
- To reduce the risk of electrical shock, regularly inspect the interconnect cables for damage.
- To reduce injury to the operator/bystander, the ClariTEE must be removed from patient contact before the application of a high-voltage defibrillation pulse.

Electromagnetic Compatibility (EMC)

Electromagnetic compatibility, or EMC, means that the device will accept interference caused by its electromagnetic environment and does not emit levels of electromagnetic energy that cause electromagnetic interference (EMI) in other devices in the vicinity. The Zura system complies with *IEC 60601-1-2 Collateral Standard: Electromagnetic Compatibility – Requirements and Tests.*

Medical electrical equipment requires special precautions regarding EMC and must be installed and operated according to these instructions. It is possible that high levels of radiated or conducted radiofrequency electromagnetic interference (EMI) from portable and mobile RF communications equipment or other strong or nearby radiofrequency sources could result in performance disruption of the ultrasound system. Evidence of disruption may include:

- · image degradation or distortion
- erratic readings
- equipment ceasing to operate
- other incorrect functioning.

If this occurs, survey the site to determine the source of disruption, and take the following actions to eliminate the source(s).

- Turn equipment in the vicinity off and on to isolate disruptive equipment.
- Relocate or reorient interfering equipment.
- Increase distance between interfering equipment and your ultrasound system.
- Manage use of frequencies close to ultrasound system frequencies.
- Remove devices that are highly susceptible to EMI.
- Lower power from internal sources within facility control (such as paging systems).
- Label devices susceptible to EMI.
- Educate clinical staff to recognize potential EMI-related problems.
- Eliminate or reduce EMI with technical solutions (such as shielding).
- Restrict use of personal communicators (cell phones, computers) in areas with devices susceptible to EMI.
- Share relevant EMI information with others, particularly when evaluating new equipment purchases which may generate EMI.
- Purchase medical devices that comply with IEC 60601-1-2 EMC Standards

Biocompatibility Safety

Biological evaluation of medical devices is performed to determine the potential toxicity resulting from contact of the component materials of the device with the body. The Zura has been designed and tested in accordance with the applicable requirements contained in ISO 10993. This is the internationally-recognized standard for the biocompatible evaluation of medical devices.

Acoustic Output

The Zura system is a Track 1 device. This means the acoustic output does not exceed the highest levels known by the Food and Drug Administration as of 1976. This is in contrast with Track 3 ultrasound diagnostic devices that contain higher acoustic output levels.

Table 10 Accepted Acoustic Output Levels for Track 1 Devices

Use	ISPTA.3 (mW/cm ²) ¹	ISPPA.3 (mW/cm ²) ²	MI ³
Cardiac	430	190	1.9

¹ISPTA.3 Derated Spatial-Peak Temporal-Average Intensity

²ISPPA.3 Derated Spatial-Peak Pulse-Average Intensity

³MI Mechanical Index

Note: Cardiac use includes transesophageal use for visualization of the heart.

Global Maximum Derated ISPTA, MI, and TI Values

ISPTA is the intensity, spatial-peak temporal-average defined as the value of the temporal-average state that ISPTA and MI are defined parameters in the FDA guidance document for ultrasound systems. TI is a temperature index defined in IEC60601-2-37, an internally recognized standard for ultrasound diagnostic equipment.

Table 11 indicates for each operating mode the values of the Thermal Index (TI), Mechanical Index (MI), derated spatial-peak temporal-average Intensity (I_{SPTA.3}), and derated peak average Intensity at peak MI (Ipa .3@MI_{max}).

Transducer Model	l _{SPTA.3} (mW/cm ²)	ТІ Туре	TI Value	мі	IPA.3 @MI _{MAX} (W/cm ²)
ClariTEE – B–6 MHz	162	Tissue	0.59	0.92	159
ClariTEE– B–7 MHz	56	Tissue	0.56	0.57	98
ClariTEE– Color Flow	63	Tissue	0.45	0.35	13

Table 11 Acoustic Output Parameters

The global maximum values for MI and TI are never greater than 1.0. As a TI or MI greater than 1.0 never occurs, neither is displayed. The user can not directly adjust the power output. The device is inherently limited to the track 1 limits and the global maximum values.

Acoustic Measurement Precision and Uncertainty

All ultrasound diagnostic devices contain a residual risk that:

- images obtained will not be clinically useful due to a lack of image quality, and
- numerical calculations such as Fractional Area Change (FAC) will be inaccurate.

These risks are mitigated by trained echocardiographers, able to discern whether an image is technically limited or clinically useful and whether quantitative assessments, such as FAC calculations, are consistent with the images. This translates into a confidence level the physician has when making a diagnosis which is reflected in his action steps. Predominate factors include the patients acoustic window quality, RFI interference, and operator technique. For the physician, the software measurement quantitative tools serve a supporting role behind the qualitative image assessment.

Labeling Symbols

Table 12 Labeling Symbols and Descriptions

\sim	Alternating current
	Protective earth (ground)
\triangle	Attention, consult Accompanying documents
0	Off (power: disconnection from the mains)
	On (power: connection to the mains)
Ò	"Off" (only for a part of EQUIPMENT)
\odot	"On" (only for a part of EQUIPMENT)
*	Type B equipment
X	Type BF equipment

Chapter 6: ImaCor Approved Devices

The ImaCor Zura has been approved for use with the following devices:

Printer Sony model UP - D897

Data Transfer Module (DTM) ImaCor-supplied drive for transferring image files, patient information, cineloops, and measurements. The drive is also used to export cineloops as DICOM files. Thumb drives are most often used, but ImaCor can provide an external hard drive to assist with large patient data transfers. **Note:** Only drives provided by ImaCor can be used. General purpose drives will not function with the ImaCor Zura system.

Chapter 7: Troubleshooting and Maintenance

Troubleshooting

If you encounter difficulty with the system refer to the troubleshooting guide below, if you are unable to resolve the issue, contact ImaCor technical support at the following numbers and addresses:

Technical support: 1-877-244-0657

Technical support email: support@imacorinc.com

Caution To reduce the risk of electric shock, do not open system enclosures.

Table 13 Troubleshooting Matrix see also the Master Message listing on page 101.

Problem	Solution	Comments
Unit does not turn on	1. Check that On/Off rocker switch is on	
	Check that power cord is connected to equipment	
	3. Check that power cord is connected to building receptacle	
	4. Check building supply	
	5. Check main fuse of machine	 Fuse holder is located integral to Power inlet. Refer to ratings on fuse holder
System computer does not boot up	 Refer to solutions for above, "Unit does not turn on" 	 Mains power is On if system fans are running. LCD and Printers appear to be powered up.
	Press pushbutton switch adjacent to ecg connector to boot computer	2. The system computer is on and booted up if the keyboard is lit
System software freezes; does not respond to user actions	Reboot computer	
System does not image	1. Toggle freeze/unfreeze command button	
	2. Reposition probe	
	3. Reboot computer	
	4. Check umbilical connection	
	5. Check probe connection to handle	
No display	1. Check LCD connections	
	2. Confirm LCD power	2. Blue LED adjacent to LCD power switch
No ECG waveform on display	1. Check ECG connection on machine	
	2. Check patient ECG lead connections	
No Print function	1. Check printer connections	1. Consists of data and power cable
	2. Confirm printer power	 LCD adjacent to printer power switch should be lit.

Maintenance

ImaCor Inc. offers a preventive maintenance plan. Contact ImaCor for further information. The following user maintenance is recommended.

Zura System Cleaning Recommendations

ImaCor recommends the following cleaning instructions for the Zura system.

Caution The ClariTEE probe is for single use only and is to be disposed after single use.

Caution To reduce the risk of electrical shock, always disconnect the power supply from the system before cleaning the system.

LCD MONITOR CLEANING

- 1. Turn off the system prior to cleaning the monitor.
- 2. Apply a small amount of one of the following recommended cleaning solutions to a soft, non-abrasive cloth:
 - Water
 - Vinegar/water solution
 - Isopropyl alcohol (i.e., Windex)
 - Petroleum benzene
- 3. Move the cloth across the display in one direction, moving from the top of the display to the bottom. Computer wipes may be used only if they specifically state they are designed for LCD screens.
 - Do not apply the cleaning solution directly to the LCD.
 - Do not scratch the LCD screen.
 - Do not use paper towels to clean the LCD screen as they may cause damage and scratches.
- 4. Never use cleaning products containing the following:
 - Acetone
 - Ethyl alcohol
 - Ethyl acid
 - Ammonia
 - Methyl chloride

OPERATOR CONSOLE CLEANING

- 1. Turn off the system prior to cleaning the operator console.
- 2. Apply a small amount of one of the following recommended cleaning solutions to a soft, non-abrasive cloth:
 - Water
 - Mild detergent (PH level at or near 7) and water solution.
 - Do not apply the cleaning solution directly to the operator console.

HANDLE AND UMBILICAL

Caution Do not immerse any part of the handle and umbilical in any cleaning solution.

- 1. Disconnect from the machine prior to cleaning.
- 2. Routinely clean the handle with special emphasis on the pads that make the connection to the probe cartridge. Any smudge or dirt on these pads could affect EEPROM communication or image quality. Use an alcohol-based cleaner; we recommend alcohol and dionized water, 70% IPA. Other types of cleaner may leave a residue that could affect performance.

Other Recommended Maintenance

POWER SUPPLY CORD

Regularly inspect the AC power cord and plug for damage.

HANDLE AND UMBILICAL

Regularly inspect the handle and umbilical for damage.

Chapter 8: Master Message Listing - Alphabetical

	Problem	Solution
Server Test Archiving server communication successful.	Informational. Zura system was able to ping archiving server.	
Configuring ECG hardware Failed to initialize ECG module		Restart the system. If the problem persists, contact ImaCor technical support.
Conflict: Last Name Patient: FRED ABDEL Select value below to resolve the conflict: © Local: ABDEL © DTM: OK	While importing or exporting patient data, the system has detected an inconsistency between the local data (that saved on the ImaCor Zura system) and the DTM.	 If the name on the main system is correct, select the Local option and click the OK button. This will append the information being transferred from the DTM to the patient record on the local system. If the name on the DTM is cor- rect, select the DTM option and click the OK button. This will change the name on the Local system to match that on the DTM
		*The local system is the one that resides within the ImaCor Zura Ultrasound Machine.
Error Critical Initialization Failure. If this error persists after reboot, contact ImaCor Tech Support. Press DK and restart the system.	 The failure of one of the following initialization steps failed: Loading configuration Loading ultrasound definitions Loading ECG definitions Loading probe settings Loading ultrasound software module Loading ECG software module Connecting to the patient data- base server 	Click the OK button and restart the system. If the error persists, contact tech- nical support.
DICOM Sync X Sync only selected patient(s) Sync all patients OK	Prompt. If at least one patient record is selected in the patient list when the Sync button is pressed, you are prompted to indi- cate whether to sync only the selected patient(s) or all patients.	Make the appropriate choice and press the OK button.
Warning Disk free space is low.	Disk space for storing patient data is limited.	Delete patients or cineloops.
Error Disk space is critically low. Operation is not allowed. OK	Disk space for storing patient data is VERY low and the current opera- tion cannot be completed.	Delete patients or cineloops.

	Problem	Solution
Exporting × Patient: {A61C3876-112B-42ED-B4DE-4A9592AD8900} File: 20120702-183322_c.dcm Exporting patient data Do not disconnect DTM. Cancel	Informational; displayed when importing data from or exporting data to a DTM. The progress bar indicates the status of the opera- tion, and the user is further advised against disconnecting the DTM.	
Server Test Failed to communicate with Archiving server. OK	The Zura system failed to ping the archiving server.	Verify that the server is online and that the Zura system is configured to communicate with the server (AET, Host, port). Also, verify that the server allows communica- tions from the Zura AE.
Server Test	The Zura system failed to ping the MWL server.	Verify that the server is online and that the Zura system is configured to communicate with the server (AET, Host, port). Also, verify that the server allows communica- tions from the Zura AE.
Failed to communicate with the patient database. Error code: -602	The system isn't communicating with the patient database.	Click the OK button and restart the system. If the problems persists, contact ImaCor technical support.
Failed to export analysis package. Error: 0x00000400	 An analysis package was requested, but it cannot be saved to the USB drive. There is insufficient space on the USB drive to write the analysis package. The USB drive was removed before the analysis package was completely written. Note: the error code shown may vary and is intended for use by an ImaCor technician. 	 If there is insufficient free space on the USB drive to accommodate the analysis package: Remove the drive from the Zura system. Insert the drive in any personal computer with a USB port. Delete unnecessary files. Files may not be deleted from the Zura system. Retry the operation If the USB drive was removed pre- maturely – i.e., before the analysis package is completely written to the drive: Reinsert the drive in the Zura system. Retry the operation.
Fried to initialize!	Review the initialization screen for the specific error.	Click the OK button and restart the system. If the problems persists, contact ImaCor technical support.

	Problem	Solution
Error Imaging is disabled. You must connect a valid probe to enable imaging.	 The Unfreeze or Acquire button was pressed while an invalid probe is connected. The probe is expired. No probe is connected. The probe is associated with a different patient. The wrong patient name may have been selected. 	 End the exam. Connect a valid probe, or- Select the correct patient name from the patient list.
Importing X Patient: {4655EA3C-CBF3-4350-967B-4DF81EBCE564} File: 20120624-081302.dcf Importing patient data Do not disconnect DTM.	Informational message. Patient data is being imported from a DTM.	No action is necessary to continue importing patient date. To halt the transfer of patient information, press the Cancel button.
Server Test	Informational. Zura system was able to ping the MWL server.	
Error No external USB drive was found.	The Analysis Package button was pressed, but no USB drive is pres- ent.	 Click the OK button. Connect a USB drive to the USB port in the Zura system. From the Configuration > Help screen, again click the Analysis Package button.
Error No external USB drive was found. OK	You are attempting to export cineloop(s) as a movie, but no USB drive is detected.	 Insert a USB thumb drive. Press the Retry button.
Warning Number of loops until maximum allowed: 5	The ImaCor system is designed to store a maximum of 100 cineloops per patient. This informational message is displayed when you near the storage limit.	No immediate action is necessary. For more information, see Deleting a Cineloop, page 62.
Warning X Patient: {A5C40655-90B4-4ECE-A12A-7A46F34DBADA} Do you want to permanently delete the selected patient data?	Warning message displayed when user attempts to delete a patient.	
Warning X Patient: (45E40555:00E44ECE: 4124-7245F340BADA) Not all onebook for this patient have been sort to the DICDM terver. Are you size you want to delete the patient and all associated onebook? Yes No	Warning message displayed when user attempts to delete a patient whose cineloops have not all been sent to the archiving server. Note: This message is only dis- played if the archiving server is active in the configuration window.	

	Problem	Solution
Printer communication failure. Check the printer status and connection.	The printer and the system aren't communicating properly.	 Check that the printer is turned on. The power switch is on the front panel. Check the data (USB) cable and power cable at the back of the printer.
Probe: ABCD-1234-TEST-0004 INVALID	The wrong patient-probe combi- nation has been selected.	Compare the ClariTEE probe's serial number label with the serial number displayed at the top if the imaging preparation screen. Most often, the selected patient is cor- rect, but the wrong probe is con- nected to the system.
Probe communication error	The probe was disconnected dur- ing the patient exam or the con- nected probe is not achieving proper contact with the handle pads	Solution A • Check for bent pins on the probe's connector • Clean handle connector pads • Reconnect probe to handle • Select patient and proceed to imaging screen Solution B
		 Connect new probe to handle Select patient and proceed to imaging screen Solution C Continue current exam; how- ever, imaging functions will be disabled
Probe Error Image: Connected probe has expired! Image: The connected probe has expired! (more than 72h have elapsed since the probe was first used) First use was: Mon Nov 12 12:02:12 2012 OK	The connected probe has been used for more than 72 hours. The message includes the date the probe was first used.	The ImaCor Zura system can be used to review the patient's recorded cineloops, but not acquire new images.
Probe identification error Image: The probe could not be identified properly. Error code: -1 Try detaching and reattaching. If the problem persists, please contact technical support. Retry Cancel	A probe is detected as new but the serial number already exists in the patient's probe database. <i>or</i> — The probe can't be identified by the system.	 Click the Cancel button to clear the error message. Detach and reattach the probe. If the problem persists, con- tact Technical support.
Probe mismatch	The patient-probe combination is invalid. Either the wrongprobe has been used or an incorrect patient record has been selected.	 Click the Cancel button to clear the error message. Press the Back button to return to the Patient Information screen and select the correct patient record. Press the Next button to advance to the imaging screen

CHAPTER 8: MASTER MESSAGE LISTING - ALPHABETICAL

	Problem	Solution
Warning Image: Control of Bitth is invalid! Please supply a valid one to continue. OK	A record has been created for a new patient, but the patient's date of birth is invalid. The patient's date of birth is entered in the fol- lowing format: mm/dd/yyyy. Examples of an invalid birthdate: • 02/17/1492 invalid year • 01/25/2079 invalid year • 07/03/194 invalid year • 99/16/1939 invalid month • 11/99/1948 invalid day	 Move the cursor to the last name field Enter the patient's last name.
Warning X Free space in megabytes(MB): 39.80 Space needed in megabytes(MB): 154.56 Not enough disk space to continue the operation. Do you want to continue anyway? Yes No	 Available disk space on the local drive is insufficient to complete an import operation, <i>or</i>- Available disk space on the DTM is insufficient to com- plete an export operation. 	 Press No to cancel the import operation. You will be able to import records to the local disk that require less than the remaining disk space. Each cineloop requires roughly 20 MB. If you are transferring data to a DTM, press No to cancel. Insert a DTM with sufficient disk space and retry the operation.
Data Transfer Module (DTM) connected Image: Connected The Export and Import functions have been enabled. To export a patient to the DTM, select the patient in the patient list and click on the Export button. OK	Informational message. The data transport module is properly con- nected, enabling patient data import and export functions.	
Warning Image: Continue is missing. Please enter to continue. OK	A record has been created for a new patient, but the patient's first name hasn't been entered.	 Move the cursor to the first name field Enter the patient's first name.
Warning The following patient record already exists: FRED ABEL Do you want to overwrite it?	The current patient already exists on the destination drive.	 Press No to return to the imaging screen. Press Yes to replace the existing patient record with the current record.
Warning Image: The Height exceeds the maximum value: 9.84 ft OK	The height entered for the patient exceeds the maximum value.9.84 ft (if Imperial)3 m (if metric)	Enter a valid height value.
Warning Image: Continue is missing! Please supply the field to continue. OK	A record has been created for a new patient, but the patient's last name hasn't been entered.	 Move the cursor to the last name field Enter the patient's last name.

A variety of actions that can produce in this remeasure. Generally speaking, this means that the operation underway was interrupted. Repeat the operation that was in grapment with the operation that was in grapment with the operation underway was interrupted. Image: the bar patter's case was added durity N An additional cineloops cannot be acquired way was interrupted. Image: the bar patter's case was added durity Image: the bar patter's case was added durity of the pattern of the pattern of the pattern of the bar pattern of the bar pattern of the bar pattern of the pattern o		Problem	Solution
Wreader An additional cineloop cannot be has been saved for this patient. An additional cineloop cannot be cineloop is overwritten. By default, the oldest existing cineloop is overwritten. By default, the oldest existing cineloop for this patient an record a new one, press the CA. OK If you wish to extern the oldest existing cineloop press the CAncel button. OK If you wish to extern the oldest existing cineloop press the CAncel button. OK If you wish to return the oldest. OK If you wish to return the oldest. If you wish to extern the acquisition is cancelled to the acquisition. If you wish to delete a cineloop to the acquisition is cancelled to the acquisition. Click the Yes button to continue the patient is cancel to actinue the p	Warning Warning Image: Control of the last patient's exam was ended abruptly! OK	A variety of actions that can pro- duce in this error message. Generally speaking, this means that the operation underway was interrupted.	Repeat the operation that was in progress when the abrupt shut- down occurred.
Warning Xi Informational. Informs user of a system shutdown and asks for confirmation. • Click the Yes button to continue system shutdown. Are you sue? Yes No • Click the Yes button to continue the patient exam. Yes No • Click the Yes button to continue the examining a patient, you return to the patient list – Warning Xi Huminse the current exam. Are you sue? • Click the wish to continue the examining a patient, you return to the patient list – Warning Xi Huminse the current exam. Are you sue? While examining a patient, you return to the patient list – Yes No • Click the yes button to continue the examining a patient, you return to the patient list –	Warning X Image: Comparison of cineloops has been reached for this patient! OK	The maximum of 100 cineloops has been saved for this patient.	 An additional cineloop cannot be acquired unless an existing cineloop is overwritten. By default, the oldest existing cineloop will be overwritten. If you wish to overwrite the oldest cineloop for this patient and record a new one, press the OK button. If you wish to retain the oldest cineloop, press the Cancel button. If you wish to retain the oldest cineloop, press the Cancel button. If you wish to delete a cineloop other than the oldest: Press Cancel to end the acquisition in progress. Delete one or more of the patient's cineloops. For more information, see Deleting a Cineloop, page 62.
Image: Construction of the second	Warning Image: Constraint of the system is about to shut down. Are you sure? Yes	Informational. Informs user of a system shutdown and asks for confirmation.	 Click the Yes button to continue with the shutdown. Click the No button to continue the patient exam.
Warning The weight entered for the patient exceeds the maximum value. Enter a valid weight value. • 1000 kg (if metric) • 1000 kg (if metric) • 2204.62 lb (if Imperial) OK Warning While examining a patient, you return to the patient information screen and select another patient in progress, press the No butto • If you wish to continue the exart in progress, press the No butto Yes No Or, while examining a patient, you return to the patient information • If you wish to end the current exam. Are you sure?	Error	 An ImaCor update package has been corrupted, or- You hasve attempted to run another executable as if it were an ImaCor update package. 	Obtain and run a valid ImaCor update package.
Warning While examining a patient, you return to the patient information screen and select another patient list – If you wish to continue the examining a patient, you return to the patient list – Yes No	Warning Warning Image: The Weight exceeds the maximum value: 1000 kg	The weight entered for the patient exceeds the maximum value. • 1000 kg (if metric) • 2204.62 lb (if Imperial)	Enter a valid weight value.
screen and create a new patient	Warning XI This will terminate the current exam. Are you sure? Yes	While examining a patient, you return to the patient information screen and select another patient name from the patient list – Or, while examining a patient, you return to the patient information screen and create a new patient	 If you wish to continue the exam in progress, press the No button. If you wish to end the current exam, press the Yes button.

CHAPTER 8: MASTER MESSAGE LISTING - ALPHABETICAL

	Problem	Solution
Warning	The ImaCor system is designed to store a maximum of 100 cineloops per patient. You have reached that maximum.	Delete one or more cineloops.
You have reached the maximum number of cineloops. Press OK to acquire over oldest loop: 20120624-080443		For more information, see Deleting a Cineloop, page 62.
OK Cancel		
Analysis Package	Informational only	
Error	The analysis package couldn't be exported to the external drive.	 Make sure the external USB is correctly inserted.
Failed to export analysis package. Error: 0x00000400	 The external USB drive was unplugged before the analysis package export was completed. 	 Make sure there is at least 5 MB free space on the external USB drive.
ОК	• The external USB drive is full.	
Error X	The analysis package must be written to an external drive, and no external drive is connected to the system.	Insert an external drive into the USB port.
Chapter 9: References

Measurements (Area and Distance)

The measurement tools provided are for use by a qualified clinician. Measurements do not translate into specific physiological parameters. The ability to obtain accurate measurements is affected by complexities of cardiac geometry, the ultrasound window presented by the patient, clarity of image, operator identification of features, operator placement of calipers, and tracing.

Sources of Measurement Errors

Measurement error can result from acoustic anomalies of the body, placement of markers by the operator and machine processing, signal acquisition, conversion, processing for display, and limitations of the monitor. The monitor pixel size is approximately 0.03 cm.

Linear data is presented in cm, rounded to the nearest 0.01 cm, representing 0.2 to 0.5% of the linear dimension of a typical LV (2 cm to 5 cm).

Area data is rounded to the nearest 0.1 cm^2 , representing 1% of typical LV end diastolic area (10 cm²), appropriate since the relative bound error in area is twice the relative error bound in length.

Accuracy

METHODOLOGY OF ASSESSING ACCURACY

Phantom side-by-side measurements (area and distance) were conducted using three predicate devices and the ImaCor Zura system. A CE tissue mimicking wire target phantom was used (ATS model 535). Areas were calculated based using the height and width measurements and based on the area formula for an ellipse (A = π hw/4). The control data used as the basis for the ImaCor Zura measurement accuracy determination were the mean height and width measurements, along with the calculated areas of three predicate systems combined.

Predicate 1 – Acuson Aspen (K934915)

Predicate 2 - Philips Sonos 7500 (K980687)

Predicate 3 - GE Vivid 7 (K051449)

The mean data entered into the below table were based on 5 measurements for the ImaCor Zura system and 15 measurements for the control group; 5 for each predicate device.

Transducer	Height (cm) (measured)	Width (cm) (measured)	Area (cm²) calculated	Comments
ImaCor Zura mean	4.02	4.39	13.83	5 measurements and area calculations
ImaCor Zura SD	0.02	0.03	0.07	
Predicate mean	4.04	4.45	14.12	Control group – 15 measurements and area calculations
Predicate SD	0.02	0.03	0.14	

Table 14 Accuracy	of Distance Measurements and Area Calculations
-------------------	--

Transducer	Height (cm) (measured)	Width (cm) (measured)	Area (cm²) calculated	Comments
ImaCor Zura accuracy	0.5%*	1.3%*	2.1%*,**	
* Assumes that the control data averages of 15 measurements by the three predicate systems, represent the actual dimensions. Does not include rounding or monitor limitations. If these limitations were included, then the relative error bounds would be approximately 2%, 3%, and 5%, respectively.				
** Need not equal 0.5% + 1.3% because errors in height and width measurement may be correlated.				

Chapter 10: Specifications

Table 15 ImaCor System Specifications

Item	Specification
System dimensions	Content
Display dimensions	17-inch LCD display
Transducers	Phased array
Imaging mode	Type B mode imaging
	Color Flow mode imaging
Application	TEE imaging. The system is optimized for the transgas- tric short axis view of the left ventricle
Measurement	Distance and area measurement capability. See "Measurements (Area and Distance)," page 108
Image storage	Data transfer module (DTM)
Cables	 Power supply cord Printer power cable Printer data cable ECG leads
Peripherals	A Sony printer is provided with the system
Probe	Single use, TEE, disposable, provided sterile
Temperature, pressure, & humidity limits	Operating limits: system • 10–40°C (50–104°F), 15–95% R.H. noncondensing • 700 to 1060 hPa (0.7 to 1.05 ATM)
	 Shipping/storage limits: system -35–65°C (-31–149°F), 15–95% R.H. noncondensing 500 to 1060 hPa (0.5 to 1.05 ATM) For storage longer than 30 days store at room temperature.
	Operating limits: Probe • 10–40°C (50–104°F), 15–95% R.H. noncondensing • 700 to 1060 hPa (0.7 to 1.05 ATM)
	 Shipping/storage limits: Probe -20-60°C (-4-140°F), 0-95% R.H. noncondensing 500 to 1060 hPa (0.5 to 1.05 ATM) For storage longer than 30 days, store at room temperature
Supply source	100–120 Vac, 7A, 50/60 Hz 200–240 Vac, 4A, 50/60 Hz
	Slide switch provided to adjust voltage to either 120 Vac or 240 Vac range
Mains fuse	3 AG Slo-Blo, 7A/250 V
	Fuse drawer style compartment located adjacent to Mains power inlet. Screwdriver needed to access fuses for replacement

Chapter 11: Glossary

Terms

ClariTEE [™]	ImaCor miniature disposable TEE probe
Cineloop	Recorded ultrasound image file
Fractional area change	Change in ventricular area at end diastole to end systole
ICU	Intensive Care Unit
ImaCor Zura [™]	ImaCor imaging system
R-wave	First upward deflection of the electrocardiogram
Real-time imaging	Imaging occurring in the present moment
TEE	Transesophageal echo where ultrasound transducer is placed in esophagus or stomach
Umbilical	Cable extending from probe handle and terminating in ZIF connector for machine connection

Acronyms

AC	Alternating current
ALARA	As low as reasonably achievable
BPM	Beats per minute
DBC	Depth brightness compensation
DTM	Data transfer module
ECG	Echocardiogram
ED	End diastolic
EMC	Electromagnetic compatibility
EMI	Electromagnetic interference
ES	End systolic
FAC	Fractional area change
FDA	Food and Drug Administration
FPS	Frames per second
ICU	Intensive care unit
ISPPA	Intensity spatial-peak pulse-average
ISPTA	Intensity spatial-peak temporal-average
LCD	Liquid crystal display
LV	Left ventricular
LVEDA	Left ventricular end-diastolic area
LVESA	Left ventricular end-systolic area

CHAPTER 11: GLOSSARY

MI	Mechanical index
RFI	Radiofrequency interference
TEE	Transesophageal echocardiography
TGC	Time gain compensation
TGSAV	Trans-gastric short-axis view
USB	Universal serial bus
ZIF	Zero insertion force