



About the ImaCor Device

ImaCor's Hemodynamic Ultrasound solution combines imaging and monitoring in a portable package for ease-of-use by high acuity care teams. The result --- lifesaving course corrections in as much as 70% of patients vs standard of care. ImaCor's miniaturized transesophageal ultrasound probe, the ClariTEE, is FDA approved to be left in place for up to 72 hours, enabling direct visualization in the ICU of cardiac filling and function over time. The patented probe has been used in over 20,000 patients with no complications. ImaCor's Hemodynamic Ultrasound provides critical care teams with pre-load and contractility data for rapid decision making, transforming care for sepsis, surgery, and cardiac support patients.

Articulating Arm

Height adjustable articulating arm with cart height reaching from 44.4" to 55.4".

Ergonomic Handle & Cord Wrap

Ergonomic handle with slender design uniquely placed for easy push-and-position.

Universal Tablet Security Bracket

Secures your tablet or iPad with lockable bracket. Locking bracket comes with 2 user keys, with or without-case sliders and supports a weight of up to 3.5lbs.

Tablet Tilt

Bracket comes with tablet tilt of 60 degrees, horizontal pan of 180 degrees and rotation around the axis of 360 degrees. Position your tablet at various angles for ease of use.

Customization

Customize your cart with Branding Options and 40+ Accessories.

Best-in-Class, Durable Materials

Anodized Aluminum, Powder-Coated Aluminum, Cyclooy Plastic, fully wipeable and approved in Operating Rooms.

Concealed Cable Management

Cable management channel allows for discreet and uninterrupted charging.

Large utility box can house power cable receptacles or other power components.

Base & Casters

4-Wheel base balances the need for stability with a small footprint. Complete with premium, quiet, non-mark and locking casters, this cart is ready to roll or stay in a locked position as you need.

Quality & Services

5 Year Limited-Warranty, ISO 9001:2015 Certified, Tip Resistant as per IEC 60601, and Remote Expert Support.

DIMENSIONS AND SPECS

- Pole Articulating Arm: 44.4" to 55.4"H
- Articulating Arm Rotation: 90 degrees
- Large Utility Box: 5"H x 8.6"L x 2.3"D
- Base Length: 17.7" to 17.9"
- Tablet Size: 6.6"-9.3"H to 9.5"-12.95"
- Tablet Thickness: 0.36" to 0.74"

WHAT'S INCLUDED

- Tryten S1 Cart: T2700
- Supply Caddy: T2557
- S-Series Handle (Silver): 551280
- Tablet Holder Slide Kit: T2563
- Custom Probe Bracket: T2405
- Triple Tab Extension Cord: 575324



CAPSAHEALTHCARE